

Gynnie

OB/GYN Stretcher



Easy to use. Multi-functional. Convenient.

Gynnie combines the durability and features of Stryker's standard stretcher with OB/GYN functionality. This multi-functional stretcher is designed for pelvic exams and can also be utilised for general use.



Pelvic exam conversion

Three easy steps to convert stretcher into a pelvic exam position with self-locking stirrups.



10 cm Ultra Comfort support surface

Heat sealed seams help prevent fluid collection.



Glideaway siderails

Swings up and out for patient comfort and allows for zero transfer gap when lowered.

Retractable fifth wheel helps improve mobility.

Permanently attached parts to maximise caregivers' time.

Optional calf supports.

Dual pedestal hydraulics to help optimise stability.

Four wheel, steel ring brakes with dual side and head activators to help secure platform and accessibility.

Gynnie

OB/GYN Stretcher



Standard features

- 10 cm (4") Ultra Comfort support surface with sealed seams
- 20 cm (8") Omni Surface casters with wheel covers
- 66 cm (26") wide patient surface
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics with constant rate descent
- Dual side-mounted foot control hydraulics
- Retractable fifth-wheel steering
- Four-wheel, steel-ring brakes with dual side and head activators
- Glideaway siderails
- Integrated storage tray with oxygen bottle holder
- Permanently-attached folding foot section
- Permanently-attached foot supports
- Pop-up push handles
- Roller bumpers
- Two IV receptacles
- One-handed foot section release lever

Specifications

Model number	1061
Overall length	205.7 cm (81")
Overall width	
Siderails up	80 cm (31.5")
Siderails down	71 cm (28")
Weight capacity	227 kg (500 lb.)
Height range (to litter top)	
High	90.8 cm (35.75")
Low	56 cm (22")
Litter positioning	
Backrest	0°– 90°
Knee gatch	0°– 90°
Trend./Reverse Trend.	± 18°
Maximum angle	87.5°
Patient surface	66 x 192 cm (26" x 75.5")
Siderails	33 x 99 cm (13" x 39")
Caster diameter	20 cm (8")
Break away point from wall	173 cm (68")

Optional features

- Hoop for drainage basin
- Patient restraints
- Serving/instrument tray
- Siderail pads
- IV poles
- Permanent folding two- or three-stage at foot end
- Removable
- Upright oxygen bottle holder

Medical

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Gynnie, Ultra Comfort, Stryker. All other trademarks are trademarks of their respective owners or holders. The products depicted are CE marked in accordance with applicable EU Regulations and Directives. This material is not intended for distribution outside the EU and EFTA.



01262005EN
DLS 03/2018
2017-16225

Copyright © 2018 Stryker
stryker.com