

# iBed<sup>®</sup> Vision Clinical Dashboard



# Stryker's Financial Guarantee

## 40% fall reduction. Guaranteed.

Through Stryker's Financial Guarantee program, we guarantee a 40% reduction in bed-related falls occurring on Stryker beds equipped with iBed wireless technology and iBed Vision\*.

Stryker's iBed Wireless Financial Guarantee: Stryker guarantees a 40% reduction in bed-related falls occurring on Stryker beds equipped with iBed Wireless Technology and iBed Vision.



\* Facility must meet eligibility requirements as outlined in the iBed Wireless Fall Reduction Contract. If a reduction of 40% of bed-related falls is not recognized after one full year (based on delivery date) of product use by the facility, Stryker will refund a portion of the license fee for iBed Vision.

# iBed<sup>®</sup> Vision Clinical Dashboard

## Personalized fall prevention

**\$34,294**

is the average cost for hospitalization cause by fall injury.<sup>2</sup>

**700,000 to 1,000,000**

patients in the United States experience a fall while being treated in the hospital.<sup>1</sup>

### Advanced fall prevention for better outcomes

We know keeping your patients safe is your #1 responsibility. As your trusted partner in prevention, we're focused on developing innovations that can help reduce falls and injuries from falls.

### Introducing iBed Vision

iBed Vision is a patient-centric clinical dashboard that provides caregivers with increased visibility to bed compliance, falls protocols and bed alarm activity. It takes near real-time data pulled directly from a patient's electronic medical records (EMR)\* and associates it with appropriate protocols to ensure a safer patient experience and an easier workflow for caregivers.

\*No PHI information is transmitted.



### **It's all about the patient**

Each patient's fall risk score is different and can change throughout his or her hospital stay. That's why iBed Vision's dynamic and personalized advanced fall prevention monitoring is so critical. Featuring configurable rules for safe bed protocol compliance, iBed Vision displays near real-time, active changes based on a patient's fall-risk – allowing you to provide optimal care as a patient's fall risk increases or decreases.

### **Smart bed alerts for personalized care**

Hospitals can set pre-determined safe bed configuration requirements that align with protocols. If the safe bed configuration isn't met (e.g., the bed exit alarm is not armed for a high fall risk patient) then notifications and alerts are wirelessly sent, independently of the nurse call cable, to iBed Vision.<sup>3</sup>

### **Wireless bed violation alerts<sup>3</sup>**

Notifications are sent from the Stryker bed when bed protocols aren't being met (e.g., low height, siderails up or down, brake set, bed exit on, and head of bed angle) independently of the nurse call cable.

### **Safe bed monitoring**


Hospital-specific safe bed protocol compliance automatically adjusts as the patient's EMR fall risk status changes.

### **Color coded status icons**

Simple design with easily recognizable bed status.

### **Customizable display**

Customize colors, icons, violation alerts and text based on hospital preference.



## Fall risk indicator

High fall risk safe bed protocol - Bed parameters

- Number of side rails in the up position
  - 3
  - 2
  - 1
- Bed exit on
- Bed exit off
- Head of bed angle at or above \_\_\_\_\_
- Low bed height
- Brake on

## Fall risk indicator

Low fall risk safe bed protocol - Bed parameters

- Number of side rails in the up position
  - 3
  - 2
  - 1
- Bed exit on
- Bed exit off
- Head of bed angle at or above \_\_\_\_\_
- Low bed height
- Brake on

FFFFFF	000000	333333	666666	999999	CCCCCC	CCCC99	9999CC	666699
660000	663300	996633	003300	003333	003399	000066	330066	660066
990000	993300	CC9900	006600	336666	0033FF	000099	660099	990066
CC0000	CC3300	FFCC00	009900	006666	0066FF	0000CC	663399	CC0099
FF0000	FF3300	FFFF00	00CC00	009999	0099FF	0000FF	9900CC	FF0099
CC3333	FF6600	FFFF33	00FF00	00CCCC	00CCFF	3366FF	9933FF	FF00FF
FF6666	FF6633	FFFF66	66FF66	66CCCC	00FFFF	3399FF	9966FF	FF66FF
FF9999	FF9966	FFFF99	99FF99	66FFCC	99FFFF	66CCFF	9999FF	FF99FF
FFCCCC	FFCC99	FFFFCC	CCFFCC	99FFCC	CCFFFF	99CCFF	CCCCFF	FFCCFF


## Choose color indicators from the above table.

- \_\_\_\_\_ High fall risk
- \_\_\_\_\_ Low fall risk
- \_\_\_\_\_ Bed Exit off
- \_\_\_\_\_ Violation notification
- \_\_\_\_\_ Patient presence

## Text Options

(Limited on character count to XX)

- Room Header
- Unit name
- Violation description (siderails, brake, generic violation)



## ibed Vision powered by iBed Wireless

Available on wireless S3 MedSurg and InTouch Critical Care beds.



## Customer Support Services

### Technical Support

Our Technical Support comprises a team of professionals available to help with your product needs. Contact via phone at 1 800 STRYKER or email at [medicaltechnicalsupport@stryker.com](mailto:medicaltechnicalsupport@stryker.com).

### Stryker's ProCare Services

Every day, you count on your medical equipment to perform at its best. With ProCare Services, our people help to ensure your equipment is ready to perform when it's needed and make it easier to get the most from your investment. When an issue arises, we promise that we'll solve it — performing repairs quickly and correctly.

ProCare isn't just a service program. ProCare is a partnership you can count on to give you one less thing to worry about, and one more reason to feel confident you're doing all you can for your clinicians, staff and patients.

All ProCare agreements provide:

- Stryker-authorized service representative
- Stryker-direct factory parts
- Two-hour callback response time
- Fixed service costs up front
- Increased uptime
- 24-72 hour equipment turnaround time\*\*

\*\*Based on the provisions of the Service Agreement and the location of the product.

### Flex Financial Program

Our financial programs provide a range of smart alternatives designed to fit your organization's needs. We offer flexibility beyond a cash purchase with payment structures that can be customized to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

### References

1. Agency for Healthcare Research and Quality. Preventing Falls in Hospitals. Retrieved December 5, 2018. <https://www.ahrq.gov/professionals/systems/hospital/fallpxtoolkit/index.html>
2. "Falls Cost U.S. Hospitals \$34 billion in Direct Medical Costs." John Hopkins Medicine. Retrieved October 18, 2019. <https://www.johnshopkinssolutions.com/newsletters/falls-cost-u-s-hospitals-30-billion-in-direct-medical-costs/>
3. Per UL1069 hospitals must still provide a wired connection for patients to send nurse call request.

3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770

**stryker.com**

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: iBed, Flex Financial, ProCare, Stryker. All other trademarks are trademarks of their respective owners or holder.