

IsoTour

Gel Support Surface




Introducing our
advanced gel support surface.

IsoTour is a gel support surface designed to help prevent all stages of pressure ulcers by addressing key risk factors of skin breakdown: pressure, shear, and moisture.

Innovative gel




TruTurn
turn assist and offload



IsoTour Pump



Low Air Loss
convertibility



IsoTour cover



IsoTour

Gel Support Surface



Standard features

- Endurance cover by Dartex cover
- Foam crib and base provide stability and lateral support
- Meets flammability standards: EN 597-1, EN 5972-2, CANADA-Method 27.7 of CAN 2-4.2 M77, ITALY UNI 9175

Additional features

- Gel mattress extender pad
- TruTurn (premium)
- Spacer fabric (premium)
- D-rings
- Anchor strap (premium)

Optional cover

- Speciality cover

Specifications

Model number	2872
Overall length	200 cm (78.75") or 214 cm (84.25")
Overall width	90.2 cm (35.5")
Thickness	24.1 cm (9.5")
Safe working load	227 kg (500 lbs)

IsoTour Pump specifications

Model number	2874
Height	43.8 cm (17.25")
Width	23.5 cm (9.25")
Depth	18.4 cm (7.25")
Weight	7.4 kg (16.3 lbs)

Warranty

- 10 years crib assembly warranty

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, IsoTour and TruTurn. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



ISOGEL - BRO - 525472EN
 ISOGEL-BRO-5_25472EN
 08/2020
 2020-25472

Copyright © 2020 Stryker
 stryker.com