

stryker

Maternity solutions



Designed for safety and simplicity

The maternity unit can be a stimulating yet demanding environment. It must not only feel comfortable for mom, but also promote mother-to-infant closeness through simple and safe products. Our maternity portfolio helps allow caregivers to focus on what matters most — providing exceptional care to mothers and their newborns.



LD304 Birthing Bed

Helping to better position your patient

Intermediate siderail position

Moves siderails back and allows mom better visibility and closeness to the bassinet and baby.

Length-adjustable footrests and calf supports

Offers nearly infinite positioning for customized comfort and slides into place with an easy-squeeze handle.

Powered Skoocher® system

Moves mom closer to the perineal edge, reducing strain on you while providing proper back support for mom.

Removable Lock-Rite® foot section

Simple, safe and exceptionally easy to handle, the Lock-Rite design provides visual confirmation that it is securely locked into place.



Nära Bassinet

Closeness when it matters most

Narrow cabinet design

A narrow cabinet design helps to promote mother-to-infant closeness from the bedside.

Two-way drawers

Offers easy access to storage, helping limit disturbance to the newborn. Soft-close slides secure drawers during transport with magnetic catch.

Height adjustable

Height adjustability provides ergonomic positioning for caregivers, and mother-to-infant closeness, regardless of bed height.

Soft-edge basket

The basket, secured to the base by three points of contact, helps provide mothers and caregivers peace of mind.



Gynnie OB/GYN Stretcher

Intuitive patient positioning and functionality

Our 4-inch UltraComfort support surface

Heat sealed seams help prevent fluid collection.

Pelvic exam conversion

Three easy steps to convert stretcher into a pelvic exam with self-locking stirrups. The permanently attached components save time and effort.

Glideaway siderails

Swing up and out for mom's comfort and provide zero transfer gap when lowered.

Retractable fifth wheel

Retractable fifth wheel helps improve mobility during transport.



Unity furniture line

Customizable pieces with a clinical edge

Our robust patient room furniture offerings marry the insights we glean from caregivers with the newest technology to create durable, cleanable, comfortable pieces for patients and users – helping to promote a healthy environment.

Bed to sleeper in seconds

Our Unity Sleeper Bench easily converts from sofa to bed with a simple pull and flip design, offering a cleanable, moisture barrier sleep surface. Six models offer the ability to customize back, armrest and tablet configurations.

Rocking motion

The Unity Glider offers a smooth, low-range rocking motion to help users relax. The curved armrests help provide comfort and support and assist users in and out of the chair safely.

Superior comfort, clinical support

Flip-up arm design on the Unity Clinical or Treatment Recliners allow caregivers to get closer to patients during patient transfers. The Performance cushion material creates a supportive and comfortable seat for both patients and users.

Customer support services

Technical support

Our technical support team comprises professionals available to help with your products from Stryker. Contact via phone at 1 800 STRYKER or email at medicaltechnicalsupport@stryker.com.

Stryker's ProCare Services

Every day, you count on your medical equipment to perform at its best. With ProCare Services, our people help to ensure your equipment is ready to perform when it's needed and make it easier to get the most from your investment. When an issue arises, we promise that we'll work to solve it — performing repairs quickly and correctly in accordance with the terms and conditions of the ProCare services agreement.

ProCare isn't just a service program. It's a partnership you can count on to give you one less thing to worry about, and one more reason to feel confident you're doing all you can for your clinicians, staff and patients.

All ProCare agreements provide:

- Stryker-authorized service representative
- Stryker-direct factory parts
- Two-hour callback response time
- Fixed service costs upfront
- Increased uptime
- 24-72 hour equipment turnaround time*

*Based on the provisions of the Service Agreement and the location of the product.

Flex Financial Program

Our financial programs provide a range of smart alternatives designed to fit your organization's needs. We offer flexibility beyond a cash purchase with payment structures that can be customized to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Cub, Flex Financial, Glideaway, Gynnie, Lock-Rite, ProCare, Nāra, Spirit Select, Stryker, Unity. All other trademarks are trademarks of their respective owners or holders.

3800 E. Centre Avenue
Portage, MI 49002 U.S.A.
t: 269 329 2100
toll free: 800 327 0770

stryker.com