

Nära®

Bassinet



When closeness matters most.

There are few moments in life more important than the time a mother first spends with her newborn.

These first moments of an infant's life are crucial to lifelong development. It is during this time that a mother bonds skin-to-skin with her newborn, creating warmth and establishing an emotional connection.

A new idea in the world of bassinets, Nära's simple, safe design brings mothers close to their newborns during the time when closeness matters most.

Key features

- Carefully designed features help allow for mother-to-infant closeness
- Simple design helps ensure ease of use and mobility for caregivers and helps to enhance interaction for mothers
- Constructed and engineered for ease of cleaning and dependable durability



Narrow cabinet design helps to promote mother-to-infant closeness from the bedside.



Height adjustable



Soft-edge basket with tilting feature



Two-way drawers

Nära Bassinet



Standard features

- Height-adjustable bassinet
- Narrow storage cabinet
- Basket secured to the base by three points of contact
- Soft-edge basket with tilt bar (7.5 ± 2.5°)
- Two-way drawers
- Sliding retractable work surface
- Open-style chart holder
- Steer locking caster
- Two caster locks
- 5" dual soft-casters
- Phthalate-free mattress cover
- Basket composed of BPA-free Triton polymer plastic
- Bassinet made with cleanable ABS plastic

Specifications

Model Number	4402-000-001 (domestic)
	4402-000-002 (international)
Overall length (handle to handle)	42.5" (108 cm)
Cabinet width	13.75" (35 cm)
Base length	33" (84 cm)
Base width	21" (53 cm)
Floor to cradle height	30.5" to 39.5" (77.5 cm to 100 cm)
Overall basket	26.5" x 12.5" x 7.5" (67 cm x 32 cm x 19 cm)
Overall mattress	26" x 12" x 1.5" (66 cm x 30.5 cm x 3.8 cm)

Stryker reserves the right to change specifications without notice.
Conversions are calculated before rounding.

Warranty

- One-year parts, labor and travel