

# Nära

## Bassinet

### When closeness matters most

There are few moments in life more important than the time a mother first spends with her newborn.

These first moments of an infant's life are special - it is during this time that a mother bonds skin-to-skin with her newborn, creating warmth and establishing an emotional connection.

A new idea in the world of bassinets, Nära's simple design brings mothers close to their newborns during the time when closeness matters most.



### Height adjustable




- Carefully designed features help allow for mother-to-infant closeness.
- Simple design helps ensure ease of use and mobility for caregivers and helps to enhance interaction for mothers.
- Constructed and engineered for ease of cleaning and dependable durability.

### Soft-edge basket with tilting feature



### Two-way drawers



# Nära Bassinet

## Standard features

- Height-adjustable bassinet
- Narrow storage cabinet
- Basket secured to the base by three points of contact
- Soft-edge basket with tilt bar ( $7.5 \pm 2.5^\circ$ )
- Two-way drawers
- Sliding retractable work surface
- Open-style chart holder
- Steer locking caster
- Two caster locks
- 12.7 cm (5") dual soft-casters
- Phthalate-free mattress cover
- Basket composed of BPA-free Triton polymer plastic
- Bassinet made with cleanable ABS plastic
- **One year warranty**



**Narrow cabinet design** helps to promote mother-to-infant closeness from the bedside.

## Specifications

Model number	4402-000-002
Overall length (handle to handle)	108.0 cm (42.5")
Cabinet width	34.9 cm (13.75")
Base length	83.8 cm (33")
Base width	53.3 cm (21")
Floor to cradle height	77 - 100 cm (30.5" - 39.5")
Overall basket	67 x 32 x 19 cm (26.5" x 12.5" x 7.5")
Overall mattress	66 x 30 x 4 cm (26" x 12" x 1.5")

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders. For more information please contact your local Stryker sales representative.

**Stryker Australia Pty Ltd**  
 8 Herbert Street  
 St Leonards NSW 2065  
 61 2 9467 1000  
[www.stryker.com.au](http://www.stryker.com.au)

**Stryker New Zealand Limited**  
 PO Box 17-136  
 Greenlane, Auckland 1546, New Zealand  
 64 9 573 1890  
[www.stryker.com](http://www.stryker.com)