

# Prime Series® stretcher



## Smart stretcher. Confident care.

Welcome to a new era of patient care with our enhanced collection of Prime Series stretchers, designed with both the patient and clinician in mind. Advanced mobility options, like Big Wheel and Zoom help reduce physical strain on care team members, while electric functionality offers patients more control over their comfort. And now, with Prime Connect, our latest smart, connected stretcher—you can help enable fall protocol compliance from anywhere in your facility.



### Advanced stretcher mobility<sup>2</sup>

Big Wheel reduces start-up push force and steering efforts, while Zoom motorized drive helps decrease manual pushing.



### Electric-powered lift pedal\*

Eliminates manual pumping, allowing for quick patient access. Simply step on the hydraulic pedal to raise the litter.



### Enhanced footboard display

The all new Prime Connect footboard includes bed exit controls, a new power and location icon, a customizable screen, a strip light to visualize bed exit status, an integrated scale locked in kilograms with the ability to view pounds and additional caregiver controls.

### Recovery chair position

Helps eliminate the need for unnecessary transfers by providing a comfortable seated position.

### Lift Assist® backrest

Uses the patient's weight to raise the head of the stretcher while the drop seat design helps reduce boosting. Also available as a fully electric option.

### IsoFlex® SE gel support surface<sup>1</sup>

Offers pressure redistribution and comfort as a component to aid in the treatment and prevention of pressure injuries.

### Power-washable\*

700 lb. safe working load.



# Prime Series stretcher

## Standard features

- Less than 21-inch low height\*
- Litter-mounted oxygen bottle holder
- Power-washable design\*
- Center-locking brake system with dual end activators
- Dual, side-mounted foot control hydraulics with uni-lower pedal
- Dual end brake/steer controls
- Glideaway siderails
- Four integrated IV receptacles
- Trendelenburg/reverse Trendelenburg positioning
- 8-inch Omni surface casters with integrated wheel covers
- Pneumatic backrest/stationary foot
- Thermoformed ABS base

## Warranty

Two years parts, labor and travel

## Electric functionality (30 inch only)

- Electric patient controls
- Foot end operator controls
- Electric lift\*
- Bed exit
- 120V or 240V available

## Prime Connect features (30 inch only)

- Nurse call integration
- Enhanced footboard
- Wireless capabilities
- Bed exit

## Optional features

- Support surfaces:
  - IsoFlex SE Support Surface
  - ComfortGel SE Support Surface
  - UltraComfort SE Support Surface
- Four-sided brake/steer controls\*
- Pop-up push handles
- Electric powered knee gatch
- Lift Assist backrest design with drop seat
- Integrated scale system
- Bed exit

## Specifications

### Model number

	26-inch litter	30-inch litter
Big Wheel advanced mobility	1115	1115
Fifth Wheel mobility	1105	1105
Zoom motorized drive system	1125	1125

### Overall length

	26-inch litter	30-inch litter
1105 and 1115	85" (215.9 cm)	85" (215.9 cm)
1125	86" (218.4 cm)	86" (218.4 cm)

### Overall width

	26-inch litter	30-inch litter
Siderails up	34" (86.4 cm)	38" (96.5 cm)
Siderails down	26" (66 cm)	30.5" (77.5 cm)
Siderails down (1125)	26" (66 cm)	31" (78.7 cm)

Patient surface	26" x 75.25" (66 x 191 cm)	30" x 75.25" (76 x 191 cm)
-----------------	----------------------------	----------------------------

Weight capacity	700 lb (318 kg)	700 lb (318 kg)
-----------------	-----------------	-----------------

### Height range

	26-inch litter	30-inch litter
High	34" (86.4 cm)	34" (86.4 cm)
Low (1105 and 1115)	20.75" (52.7 cm)	20.75" (52.7 cm)
Low (1125)	23" (58 cm)	23" (58 cm)

### Litter positioning

	26-inch litter	30-inch litter
Backrest	0°– 90°	0°– 90°
Knee gatch	0-5.5"	0-5.5"
Trend./reverse Trend. (1105 and 1115)	+17°/-17° (± 3°)	+17°/-17° (± 3°)
Trend./reverse Trend. (1125)	+17°/-16° (± 3°)	+17°/-16° (± 3°)

Siderail length	14" x 57.75" (36 x 147 cm)	14" x 57.75" (36 x 147 cm)
-----------------	----------------------------	----------------------------

Caster diameter	8" (20 cm)	8" (20 cm)
-----------------	------------	------------

EMC	IEC 60601-1-2	IEC 60601-1-2
-----	---------------	---------------

- Three-sided hydraulic controls\*
- Dual-ended siderail release
- Permanent two-stage IV pole
- Serving tray
- Defibrillator tray
- Defibrillator tray/foot extender/charting surface
- Chartholder
- Serving tray holder/footboard
- Siderail pads
- IV caddy
- Upright oxygen bottle holder
- Restraint straps

\* 1105 and 1115 models only

## References

1. Recommended to be implemented in combination with clinical evaluation of risk factors and skin assessments made by health care professional.
2. Stryker internal data report. 2017.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Big Wheel, ComfortGel, Glideaway, IsoFlex, Lift Assist, Prime Series, Stryker, Zoom. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.