

stryker

ProCuity™ bed series



Brilliance in a bed.

Our most intelligent bed yet, ProCuity features low height for all patient care environments, advanced fall prevention technologies¹, cable-less nurse call and is a truly wireless bed platform*. ProCuity offers a proactive approach to help prevent falls, reduce injury from falls, and streamline caregiver workflow¹.



L model

29cm low height
10cm touchscreen



LE model

29cm low height
Electric brakes
USB outlet



Z model

36cm low height
Zoom Motorised Drive



ZM model

Premium 20cm screen
Powered mattress
integration ability



Working across departments as one seamless support system, ProCuity offers caregivers one single platform – helping to make standardisation across your hospital easier and improve caregiver efficiencies.

• iBed Wireless and Cable-less nurse-call ready •

Specifications

ProCuity Specifications

Model number 3009

Overall length 228.6 cm

Width (siderails up)
105.4 cm

Width (siderails down)
100.3 cm

Height range (to litter top) L & LE models

High: 76.2 cm

Low: 29.2 cm

Height range (to litter top) Zoom versions

High: 81.3 cm

Low: 35.6 cm

Litter positioning

Backrest: 0° to 65° ±5°

Knee gatch: 0° to 30° ±5°

Trend./Reverse Trend.: +12° to -12° ±5°

Patient Surface

213.4 cm x 88.9 cm)

Patient Surface with bed extender:

243.8 cm x 88.9 cm)

Caster Diameter: 15.2 cm (6")

Max. Patient Weight: 226.8 kg²

1. Design Powered by Research: Improving Patient and Caregiver Safety Through Evidence-Based Design. Guy Fragala. August 2020. 2. Stryker data held on file.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders. For more information please contact your local Stryker sales representative.

Standard Features

- Bed scale**
- 3-zone Adaptive Bed Alarm
- iBed Watch
- Auxiliary outlets 230 VAC, 50 Hz, 8A
- Battery back-up
- Integrated bed extender
- Nurse call interface
- Pendant port
- Manual brakes
- Patient belongings holder
- Wireless capability*
- Motion interrupt
- Integrated Foley hook
- Integrated restraint hooks
- Automatic under-bed light
- Integrated Pump holder
- 3-position siderails
- Smart Prompts and Safety Dashboard

Electronic Specs

- Agency approvals: UL approved to IEC 60601
- Volts: 100-240 VAC
- Ampere rating: 8A
- Current leakage: Not more than 300 μ A (microamps)
- Frequency: 50/60 Hz

*only available on the LEX, ZX & ZMX models. No patient information is transferred. **1.4 kg (\pm 3 lb) of total patient weight for patients who weigh 27.2 kg (60 lb) to 45.4 kg (100 lb)