

Symmetry® II

Recliner Bed



Easily converts to lounge or sleep positions

High-resiliency polyurethane foam and fiber-wrap cushioning provide maximum comfort and lasting shape

Equipped with front-locking casters

Limited lifetime warranty

Premium Strength, Quality Care

With Stryker's Recliner Bed, expect durability without compromising comfort. Three highly-versatile positions (chair, lounge and full flat-bed) will help guests feel right at home in any patient room.



Lounge Position



Full Flat-bed Position

Symmetry® II

Recliner Bed

stryker®



Medical

3800 E. Centre Ave.
Portage, MI 49002
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- 3" (8 cm) back swivel casters
- 3" (8 cm) fixed front locking casters
- CAL 117

Specifications

• Model	Recliner Bed	
	Part #3500-000-650	
• Overall Width	35"	(89 cm)
• Overall Height	35.5"	(90 cm)
• Depth		
Overall	37"	(94 cm)
Seat	21"	(53 cm)
Lounge Position	66"	(168 cm)
Flat-bed Position	85"	(216 cm)
• Weight Capacity	300 lbs	(136 kg)
• Weight	130 lbs	(59 kg)
• Seat Width	27"	(69 cm)
• Seat Height From Floor	17.5"	(45 cm)
• Arm Height From Floor	26"	(66 cm)
• Overall Sleep Area	27" x 78"	(69 cm x 198 cm)

*Stryker reserves the right to change specifications without notice.
Furniture is manufactured exclusively for Stryker Medical.*

Warranty

- Limited lifetime warranty

Optional Features

- Choose from a variety of fabrics and stains
- CAL 133

Copyright © 2006 Stryker
PC/SCG-7.5M-0508, 06-SMD-1292
Mkt Lit-94-02022005 Rev D

