

Spirit Plus[®] Bed



A low-height solution to meet your essential needs

When paired with strategic hospital processes, the Spirit Plus Bed can help create a safer hospital environment for both patients and staff.

The implementation of Spirit Plus beds can help reduce per patient costs and improve patient outcomes. Boasting a low bed height of 10.75 inches, a night light and siderails designed to provide optimal support, the Spirit Plus bed helps keep patients safe when paired with hospital protocols. Quad siderail design with specialty short head rails can help support throughout patient transfers, helping reduce the risk to their safety.



Siderail-embedded controls allow patients and staff to easily adjust head and knee articulations.



Open architecture litter deck compatible with most standard mattresses.



Quad siderail design with specialty short head rails can help support throughout patient transfers reducing risk to their own safety.

10.75-inch low bed height allows the patient to enter and exit with stability.

Open architecture allows facilities to customize the mattress used to the needs of each patient.

Integrated nurse call system provides better staff and patient communication.

Central lock-and-steer brake function helps increase caregiver control.

Spirit Plus Bed



Standard features

- Battery back-up
- Patient lock-outs
- Dual lock-and-steer brake system
- Footboard control panel
- Auto contour
- Low chair positioning
- Trendelenburg/reverse Trendelenburg
- Emergency CPR — automatic
- Total bed lock-out
- Molded plastic headboard/footboard
- Embedded patient controls in siderails
- Integrated nurse call button
- Power status indicator
- High quad rails

Loaded model

Includes all features listed above, plus:

- 120 V AC accessory outlet
- Obstruction sensors
- 4-inch tool-less bed extension

Specifications

Model number	5600
Overall length	91.75" (233 cm)
Overall width	40" (101.6 cm)
Safe working load	500 lb. (227 kg)
Height range (to litter top)	
High	34" (86 cm)
Low	10" (27.3 cm)
Litter positioning	
Knee gatch	0°–25°
Trend./Reverse Trend.	14° (+/-)
Backrest	0°–60°
Foot section	0°–10°
Patient surface	
Minimum	35" x 80" (89 x 203 cm)
Maximum	35" x 84" (89 x 213 cm)
Caster diameter	5" (12.7 cm)

Electronics

Volts	120 VAC
--------------	---------

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Spirit Plus, Stryker. All other trademarks are trademarks of their respective owners or holder.