


Spirit[®] Behavioral Health Platform Bed



Contemporary design with robust construction

Our behavioral health platform bed combines contemporary design with robust construction. Tamper-resistant hardware, durable frame materials and a large weight capacity make this bed ideal for behavioral health facilities requiring a resilient bed frame. The flat surfaces of the bed make it easy to clean and sanitize.

Headboard and footboard are high density polyethylene (HDPE)

- Easy to clean; stain and odor resistant
- Will not chip, swell, rot or crack

Tamper-resistant fasteners used throughout bed

Standard model (80")

- Restraint loops (optional)
- Rubber feet
- 3-inch headboard / 3-inch footboard
- Bottom pan

Extended model (84")

- Restraint loops
- Rubber feet
- 7-inch headboard / 3-inch footboard
- Bottom pan

Specifications	Standard model	Extended model
Model number	5800	5800
Overall length	81.5" (207 cm)	85.5" (217.2 cm)
Mattress length	80" (203.2 cm)	84" (218.3 cm)
Overall width	36" (91.4 cm)	36" (91.4 cm)
Weight capacity	1000 lb. (454 kg)	1000 lb. (454 kg)
Height (to liter top)	16" (40.6 cm)	16" (40.6 cm)

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Spirit, Stryker. All other trademarks are trademarks of their respective owners or holder.