

Spirit[®] Behavioral Health Bed



Low height, high security.

Behavioral health departments can place unique demands on hospital employees. We worked to create a solution with features that help to satisfy these demands and increase security within these specific units when partnered with specialized hospital processes. In addition to the specialized security features that can be added to Spirit beds, all beds have a 10.75-inch low bed height. This low height helps to reduce fall injuries and increase patient comfort to create a more secure sleeping environment.



Any or all patient controls can be locked to restrict use with the patient lockout buttons. Behavioral health beds are equipped with an additional security code lock-out.



Optional bed exit alarm monitors patient's bed exit activity in single-zone or three zone.



Nonremoveable mattress deck.

Siderails comply with the FDA's monitored Entrapment Zones 1 – 4

Filled siderails, headboard and footboard to help prevent harm

Patient lockout restricts all bed functions

Nonremoveable mattress deck, headboard, and footboard

Tamper-resistant hardware

Permanent control box cover

Removable brake pedal

Shortened power cord

Spirit Behavioral Health Bed



Spirit Plus Behavioral Health Bed

- Battery back-up
- Patient lock-outs
- Dual lock-and-steer brake system
- Footboard control panel
- Auto contour
- Low chair positioning
- Trendelenburg/reverse Trendelenburg
- Total bed lock-out
- Molded plastic headboard/footboard
- Embedded patient controls in siderails
- Integrated nurse call button
- Power status indicator
- High quad rails
- Filled rails and boards
- 4-inch tool-less bed extension
- Behavioral health package (tamper-resistant hardware, nonremovable headboard and footboard, removable brake pedal)
- Security lockout
- Shortened power cord

Spirit Select Base Behavioral Health Bed

includes all features listed above, plus:

- Backlit staff controls
- Level-all bed reset
- Built-in scale
- Integrated single-zone bed exit alarm (optional 3-zone alarm)
- Programmable night light

Spirit Select Loaded Behavioral Health Bed

includes all features listed above, plus:

- Watchdog Bed Status Indicator
- Obstruction sensors
- Nurse reminder
- Additional colored night lights
- 120V AC accessory outlet

Specifications

Model number	5600/5700 (plus/select)
Overall length	91.75" (233 cm)
Overall width	40" (101.6 cm)
Safe working load	500 lb. (227 kg)
Height range (to litter top)	
High	34" (86 cm)
Low	10.75" (27 cm)
Litter positioning	
Knee gatch	0°–25°
Trend./Reverse Trend.	14° (+/-)
Backrest	0°–60°
Foot section	0°–10°
Patient surface	
Minimum	35" x 80" (89 x 203 cm)
Maximum	35" x 84" (89 x 213 cm)
Caster diameter	5" (12.7 cm)

Electronics

Volts	120 VAC
--------------	---------

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

stryker.com

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Spirit, Stryker, Watchdog. All other trademarks are trademarks of their respective owners or holder.