


Spirit Select[®] VA Bed



Helping to create a home-like environment in low height.

Spirit Select beds are used widely throughout the VA healthcare system. Patient falls and safe patient handling are two important VA initiatives, and, in concurrence with proper hospital protocols, the Spirit Select has helped many VA hospitals achieve their safety goals. With a low height of 10.75 inches, Spirit Select can help reduce patient falls and injuries.

This low height also makes it easier for patients to enter and exit bed, can help reduce falls during transfers, and maximizes staff transfer assistance. Best of all, the Spirit Select's Bed Exit Alarm protects staff and patients from false alarms, alarm desensitization and nonrestoration of settings.



10.75 inch low bed height can help patients enter and exit the bed with stability.



Watchdog Bed Status Indicator continuously monitors distinct bed data points.



Conveniently positioned buttons reduce staff back strain and minimizes lifting when caring for patients.

Bed Exit Alarm fatigue prevention helps protect staff and patients from common alarm hazards.

Customizable bumper lights are programmed to turn on when patient exits the bed to help reduce trips and falls.

Built-in scale reduces need to transfer patient for weighing.

Spirit Select VA Bed



Standard features

- Battery back-up
- Patient lock-outs
- Dual lock-and-steer brake system
- Footboard control panel
- Auto contour
- Low chair positioning
- Trendelenburg/Reverse Trendelenburg
- Emergency CPR — automatic
- Emergency CPR — manual
- Total bed lock-out
- Molded plastic headboard/footboard
- Embedded patient controls in siderails
- Integrated nurse call button
- Power status indicator
- Level-all bed reset
- 4-inch tool-less bed extension (option)
- 120 VAC accessory outlet
- Obstruction sensors (option)
- High quad rails
- Watchdog Bed Status Indicator
- Wood headboard and footboard
- Woodgrain inlays

Specifications

Model number	5700
Overall length	91.75" (233 cm)
Overall width	40" (101.6 cm)
Safe working load	500 lb. (227 kg)
Height range (to litter top)	
High	34" (86 cm)
Low	10.75" (27 cm)
Litter positioning	
Knee gatch	0°–25°
Trend./Reverse Trend.	14° (+/-)
Backrest	0°–60°
Foot section	0°–10°
Patient surface	
Minimum	35" x 80" (89 x 203 cm)
Maximum	35" x 84" (89 x 213 cm)
Caster diameter	5" (12.7 cm)

Electronics

Volts	120 VAC
--------------	---------

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Spirit Select, Stryker, Watchdog. All other trademarks are trademarks of their respective owners or holder.

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

patientcare.stryker.com

Mkt Lit-1053 08 JAN 2015 Rev C.2
Copyright © 2017 Stryker
Printed in U.S.A