

Spirit Select[®] Youth Inlay Bed




Helping to create a safe environment for children.

The Spirit Select Youth Inlay Bed helps to ease the transition between home and the hospital environment by allowing young patients to rest more comfortably during their stay.

The low-height, innovative quad siderails and a programmable night light have all been designed to help patients enter and exit bed properly and safely.

Open architecture, easy-to-clean bed frame and adjustable positions can help caregivers tend to their patients' needs while reducing risk to their own safety.



10.75 inch low bed height helps reduce fall injuries and encourages children to enter and exit bed properly and safely.



Patient controls can be locked to restrict use with the patient lock-out buttons.



Playground inlay graphics helps create a feeling of comfort for the youth patients and their families.

Customizable bumper lights with the adjustable night light can help children feel more comfortable at night.

Obstruction sensors helps prevent bed from lowering and causing injury or damage.

Single-zone bed exit alarm monitors a child's bed exit activity

Spirit Select Youth Inlay Bed



Standard features

- Battery back-up
- Patient lock-outs
- Dual lock-and-steer brake aystem
- Footboard control panel
- Auto contour
- Low chair positioning
- Trendelenburg/Reverse Trendelenburg
- Emergency CPR — automatic
- Emergency CPR — manual
- Total bed lock-out
- Molded plastic headboard/footboard
- Embedded patient controls in siderails
- Integrated nurse call button
- Backlit staff controls
- Power status indicator
- Level-all bed reset
- Built-in scale
- Integrated single-zone bed exit alarm
- Programmable night light
- High quad rails
- 4-inch tool-less bed extension (option)
- 120 VAC accessory outlet
- Obstruction sensors (option)
- Youth siderail inlays (playground theme)

Specifications

Model number	5700
Overall length	91.75" (233 cm)
Overall width	40" (101.6 cm)
Safe working load	500 lb. (227 kg)
Height range (to litter top)	
High	34" (86 cm)
Low	10.75" (27 cm)
Litter positioning	
Knee gatch	0°–25°
Trend./Reverse Trend.	14° (+/-)
Backrest	0°–60°
Foot section	0°–10°
Patient surface	
Minimum	35" x 80" (89 x 203 cm)
Maximum	35" x 84" (89 x 213 cm)
Caster diameter	5" (12.7 cm)

Electronics

Volts	120 VAC
--------------	---------

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Spirit Select, Stryker. All other trademarks are trademarks of their respective owners or holder.

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

patientcare.stryker.com

Mkt Lit-1052 08 JAN 2015 Rev C.2
Copyright © 2017 Stryker
Printed in U.S.A