

# Stretcher Chair

## Versatile. Simple. Compact.

Designed for surgery centres and small spaces, Stryker's stretcher chair benefits patient throughout without compromising patient care. Choose positioning from horizontal to upright to achieve desired performance and patient comfort.



### Eye Stretcher Chair

The eye stretcher chair provides all the benefits of the standard stretcher chair plus a dual-articulation, enhanced clearance headpiece for surgical positioning.

Stretcher easily articulates from **flat to a chair position**.

Pneumatic backrest manoeuvres **0 to 90 degrees**.

**Dual hydraulic jacks and four wheel, steel-ring brakes** enhance stability.

**Foot section can be moved dependently or independently** of backrest.

**Multiple sticker colour options** allow you to customise the stretcher for identification purposes.



Grey



Red



Blue

# Stretcher Chair



## Standard features

- 5 cm (2") vinyl mattress
- 15 cm (6") casters
- 61 cm (24") wide patient surface
- Adjustable height footrest
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual-side activators
- Hands-free Trendelenburg
- Swing-down siderails with comfort armrests
- Thermoformed ABS hood with integrated storage tray and oxygen bottle holder
- Multi-position pushbar (model 5050 only)
- Independent foot section
- Dual-articulation, enhanced-clearance headpiece (5051)
  - 8 cm (3"), 5 cm (2") or 2.5 cm (1") flat

## Specifications

<b>Model number</b>	5050/5051
<b>Overall length</b>	193 cm (76")
<b>Width</b>	76 cm (30")
Overall	76 cm (30")
Litter	70 cm (27.5")
<b>Weight capacity</b>	181 kg (400 lb.)
<b>Height range (to litter top)*</b>	
High	85 cm (33.5")
Low	55.5 cm (22")
<b>Litter positioning</b>	
Backrest	0° - 90°
High	0° - 80°
Low	± 18°
<b>Siderails</b>	± 45°
Standard	25 x 79 cm (10" x 31")
<b>Patient surface</b>	61 x 188 cm (24" x 74")
<b>Caster diameter</b>	15 cm (6")

## Optional features

- 5 cm (2") Enhanced Comfort Mattress
- 8 cm (3") Enhanced Comfort Mattress
- 8 cm (3") or 5 cm (2") Concave head pad (model 5051 only)
- IV poles
  - Tethered
  - Standard removable
- Restraint straps
- Extended warranty available
- Adjustable arm board and clamp (5051) (surgery accessory rail required)
- Air delivery system/drape support (5051)
- Surgery accessory rail (5051) (not available with patient transfer siderails)
- Temporal or superior wristrest (5051)

## Medical

**This document is intended solely for the use of healthcare professionals.**

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



16226EN  
DLS 02/2018  
2017-16226

Copyright © 2018 Stryker  
stryker.com