

Stretcher Chair

Versatile. Simple. Compact.

Ideal for surgery centers and small spaces, our stretcher chair maximizes patient throughput without compromising patient care. Choose optimal positioning from horizontal to upright to achieve exceptional performance and patient comfort.



Eye Stretcher Chair

The eye stretcher chair provides all the benefits of the standard stretcher chair plus a dual-articulation, enhanced clearance headpiece for surgical positioning.

Stretcher easily articulates from **flat to a chair position**

Pneumatic backrest maneuvers **0 to 90 degrees**

Dual hydraulic jacks and four wheel, steel-ring brakes optimize stability

Foot section can be moved dependently or independently of backrest

Multiple sticker color options allow you to customize the stretcher for identification purposes.



Gray

Red

Blue

Stretcher Chair



Standard features

- 2" (5 cm) vinyl mattress
- 6" (15 cm) casters
- 24" (61 cm) wide patient surface
- Adjustable height footrest
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual-side activators
- Hands-free Trendelenburg
- Swing-down siderails with comfort armrests
- Thermoformed ABS hood with integrated storage tray and oxygen bottle holder
- Multi-position pushbar (model 5050 only)
- Independent foot section
- Dual-articulation, enhanced-clearance headpiece (5051)
 - 3" (8 cm), 2" (5 cm) or 1" (2.5 cm) flat
 - 3" (8 cm) or 2" (5 cm) concave

Specifications

Model number	5050/5051
Overall length	76" (193 cm)
Width	30" (76 cm)
Overall	30" (76 cm)
Litter	27.5" (70 cm)
Weight capacity	400 lb. (181 kg)
Height range (to litter top)*	
High	33.5" (85 cm)
Low	22" (55.5 cm)
Litter positioning	
Backrest	0° - 90°
High	0° - 80°
Low	±18°
Siderails	±45°
Standard	10" x 31" (25 x 79 cm)
Patient surface	24" x 74" (61 x 188 cm)
Caster diameter	6" (15 cm)

Optional features

- 2" (5 cm) Enhanced Comfort Mattress
- 3" (8 cm) Enhanced Comfort Mattress
- Concave head pad (model 5051 only)
- IV poles
 - Tethered
 - Standard removable
- Restraint straps
- Extended warranty available
- Adjustable arm board and clamp (5051) (surgery accessory rail required)
- Air delivery system/drape support (5051)
- Surgery accessory rail (5051) (not available with patient transfer siderails)
- Temporal or superior wristrest (5051)

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Stryker**. All other trademarks are trademarks of their respective owners or holder.

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

patienthandling.stryker.com

Copyright © 2017 Stryker
Mkt Lit-64 01262005 Rev G.2
Printed in U.S.A.