



## Stretcher Chair

Ideal for surgery centres and small spaces, our Stretcher Chairs maximise patient throughput without compromising patient care and outcomes. Choose optimal positioning from horizontal to upright to enhance functionality and patient comfort.

**Pneumatic backrest** manoeuvres 0 to 90 degrees so it easily articulates from a flat to a chair position.

**Dual hydraulic jacks and fourwheel**, steel-ring brakes optimise stability.

**Foot section** can be moved dependently or independently of backrest.

**Multiple bumper colour options** allow you to customise the stretcher for identification purposes.



Navy Gray Red

## Eye Stretcher Chair



The Eye Stretcher Chair provides all the benefits of the standard Stretcher Chair plus a dual-articulation, enhanced clearance headpiece for surgical positioning.

# Stretcher Chair

Versatile. Simple. Compact.

## Standard features

- 5.1 cm (2") vinyl mattress
- 15.2 cm (6") casters
- 61.0 cm (24") wide patient surface
- Adjustable height footrest
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual-side activators
- Hands-free trendelenburg
- Swing-down siderails with comfort armrests
- Thermofomed ABS hood with integrated storage tray and oxygen bottle holder
- Multi-position pushbar (model 5050 only)
- Independent foot section
- **Two year warranty**

## 5051 Standard features (All standard features plus)

- Dual-articulation, enhanced-clearance headpiece
- Head pads - Ultra Comfort: 7.6 cm (3"), 5.1 cm (2") or 2.5 cm (1") flat; 7.6 cm (3"), 5.1 cm (2") concave

## Specifications

<b>Model number</b>	5050/5051
<b>Overall length</b>	193.0 cm (76")
<b>Width</b>	
Overall	76.2 cm (30")
Litter	69.9 cm (27.5")
<b>Weight capacity</b>	181.4 kg (400 lb.)
<b>Height range (to litter top)</b>	
High	85.1 cm (33.5")
Low	55.9 cm (22")
<b>Litter positioning</b>	
Backrest	0° - 90°
Leg	0° - 80°
Trend./reverse Trend.	±18°
<b>Patient surface</b>	61 x 188 cm (24" x 74")
<b>Siderails</b>	±45°
Standard	25 x 79 cm (10" x 31")
<b>Caster diameter</b>	15.2 cm (6")

## Optional features

- 5.1 cm (2") Enhanced Comfort Mattress
- 7.6 cm (3") Enhanced Comfort Mattress
- Concave head pad (model 5051 only)
- IV poles: tethered; standard removable
- Restraint straps

## 5051 Optional features (All optional features plus)

- Adjustable arm board and clamp (surgery accessory rail required)
- Air delivery system/drape support
- Surgery accessory rail (not available with patient transfer siderails)
- Temporal or superior wristrest

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

For more information please contact your local sales representative.

**Stryker Australia Pty Ltd**  
8 Herbert Street  
St Leonards NSW 2065  
61 2 9467 1000  
www.stryker.com.au

**Stryker New Zealand Limited**  
PO Box 17-136  
Greenlane, Auckland, New Zealand  
64 9 573 1890  
www.stryker.com

JMEH\_1974\_1\_AOT  
Copyright © 2020 Stryker