

Eye Surgery Stretcher and Surgistool® Chair



Eye Surgery Stretcher

Access. Simple. Functionality.

Stryker's Eye Surgery Stretcher offers enhanced headend clearance*, giving surgeons ample leg room to manoeuvre while performing delicate eye procedures. The reliable work surface and mobile design combine to help with surgical throughput.



Enhanced head end clearance maximises the surgeon's access for either a temporal or a superior approach to surgery.



Optional contoured or flathead piece comfortably supports patient during procedures.

Retractable fifth wheel helps improve mobility.

Metal tooth and groove braking system prevents swivel and roll.

Option of two- or three-position fold-down siderails.

Multiple bumper color options allow you to customise the stretcher for identification purposes.



Navy Gray Red

*Compared to other models.

Surgistool Chair

The ideal surgeon's stool

Surgistool Chair offers customisable seat configurations with fully adjustable backrest and armrests provide optimal support to meet the surgeon's specific positioning preferences. Easy pedal access and a stable yet mobile design combine to deliver a highly functional addition to your operating room.



Three spindle base provides superior control and pedal access.



Optional padded armrests allows versatile positioning and stability.

Multiple seat configurations offer personalised comfort.



Teardrop seat



Wedge seat



Saddle seat
(optional)

Eye Surgery Stretcher and Surgistool Chair



Standard features: Eye Surgery Stretcher

- 7.6 cm Enhanced Comfort Mattress
- 20.3 cm Omni Surface casters with wheel covers
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side hydraulic controls
- Fold-down siderails
- Four IV receptacles
- Central-locking brakes with foot-end brake and steer pedal
- Hands-free Trendelenburg
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Rubber bumpers
- Stationary foot-end litter
- Uni-lower pedal
- Dual-articulation, enhanced clearance headpiece

Standard features: Surgistool Chair

- Adjustable contoured backrest
- Dual-density foam seat cushions (except saddle seat option)
- Foot-controlled, hydraulic, three-spindle base
- Front-end locking system
- Thermoformed plastic hood
- Teardrop or wedge seat

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

For more information please contact your local sales representative.

Stryker Australia Pty Ltd
8 Herbert Street
St Leonards NSW 2065
61 2 9467 1000
www.stryker.com.au

Stryker New Zealand Limited
PO Box 17-136
Greenlane, Auckland, New Zealand
64 9 573 1890
www.stryker.com

JMEH_2078_1_AQT
Copyright © 2020 Stryker

Specifications

	Eye Surgery Stretcher	Surgistool Chair
Model number	1089	830
Overall length	228.6 cm (90")	
Overall width	80.0 cm (31.5")	
Weight capacity	226.8 kg (500 lb.)	136.1 kg (300 lb.)
Height range		
High	86.4 cm (34") (to litter top)	72.4 cm (28.5") 77.5 cm (30.5") with riser kit
Low	56.5 cm (22.25") (to litter top)	53.3 cm (21") 58.4 cm (23") with riser kit
Seat cushion style		Teardrop / wedge / saddle
Litter positioning		
Backrest	0° - 90°	
Knee gatch	0° - 30°	
Trend./Reverse Trend.	± 18°	
Patient surface	66 x 221 cm (26" x 87")	
Siderails	33 x 139.5 cm (13" x 55")	
Caster diameter	20.3 cm (8")	

Stryker reserves the right to change specifications without notice.

Optional features: Eye Surgery Stretcher

- 10.2 cm Ultra Comfort Mattress
- Pre/post-op shoulder pads assembly (crank backrest only)
- Air delivery system/drape support
- Crank backrest
- Crank knee gatch
- Defibrillator tray/foot extender
- Temporal or superior wristrest
- Footboard/chartholder
- Three-sided hydraulic controls
- Three-sided brake and steer pedals
- Integrated patient transfer/arm board
- IV caddy
- Restraint straps
- Serving/instrument tray
- Siderail pads
- Surgery accessory rail
- Three-position siderails
- Upright oxygen bottle holder
- Oxygen bottle retaining ring
- Serving/instrument tray holder
- Head pads - Ultra Comfort
 - 7.6 cm, flat
 - 7.6 cm, concave
- IV poles
 - Permanent folding two- or three-stage at foot end
 - Removable
- Extended warranty available

Optional features: Surgistool Chair

- 5.1 cm riser kit
- Saddle seat
- Multi-position armrest