

Symmetry Plus Patient Room Recliners

Stylish. Versatile.

Versatile recliners for a healing environment.

Our stylish patient room recliners maximise patient and guest comfort, while providing uncompromised durability. Choose from Rocker or Wall-Saver models, with open arm or upholstered arm options.



Backrest stops in infinite number of recline positions.

Posture-correct lumbar and lateral cushioning.

Pull-out foot tray gives patients added comfort while in transit.

Rocker model's rear lock lever allows chair to remain stationary.

Wall-saver model reclines away from the wall to save space.

Sinuuous wire-spring construction provides comfort over time.

Rigid push bar for easy manoeuvrability.

Wide arm caps are extended to facilitate easy entry and exit.

Rocker recliner

- Rear lock lever allows chair to remain stationary


Wall-saver

- Reclines away from the wall to save space

Step-lock brakes



Rear push bar



Rolling footrest



Symmetry Plus Patient Room Recliners

Stylish. Versatile.

Standard features

- Stationary or mobile
- Mobile with 7.5 cm front rigid/rear swivel locking casters
- Mid-ottoman cushion provides extra support when feet are elevated
- Chair construction includes metal seat frame and sinuous wire springs
- Removable back cushion enables thorough cleaning and service
- Highly-resilient polyurethane foam cushions
- Stationary arms and fixed sides provide stability for drainage catheters or electromechanical mounts (Open-arm only)
- Egress assist polyurethane arm caps
- Meets flammability standards of BIFMA F1-1978 – revised 1980 and CAL 117
- Meets ANSI/BIFMA X5.4 2005 per tests conducted in a certified test facility
- Rigid push bar for easy manoeuvrability
- Pull-out rolling footrest gives patient added comfort while in transit
- **Limited lifetime warranty**

Optional features

- Open arm or upholstered arm
- Egress assist wood arm caps
- Hinged push bar for manoeuvrability
- Four-wheel swivel with locking casters
- Choose from a variety of stains and fabrics
- CAL 133 in selected fabrics
- Foley Bag hooks (Wall-saver only)

Specifications

Model	Rocker Recliner 3500-000-750	Wall-saver 3500-000-730
Overall width	78.7 cm (31")	78.7 cm (31")
Overall height	115.6 cm (45.5")	120.7 cm (47.5")
Depth		
Overall	105.4 cm (41.5")	102.9 cm (40.5")
Seat	49.5 cm (19.5")	49.5 cm (19.5")
Recline/legrest up	177.8 cm (70")	177.8 cm (70")
Weight	62.6 kg (138 lb.)	59.2 kg (130.5 lb.)
Weight capacity	158.8 kg (350 lb.)	158.8 kg (350 lb.)
Seat width	58.4 cm (23")	58.4 cm (23")
Back height from seat	82.6 cm (32.5")	85.1 cm (33.5")
Seat height from floor	52.1 cm (20.5")	52.1 cm (20.5")
Arm height from floor	68.6 cm (27")	68.6 cm (27")

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

For more information please contact your local sales representative.

Stryker Australia Pty Ltd
8 Herbert Street
St Leonards NSW 2065
61 2 9467 1000
www.stryker.com.au

Stryker New Zealand Limited
PO Box 17-136
Greenlane, Auckland, New Zealand
64 9 573 1890
www.stryker.com