### stryker

0

# Transport Stretcher

### **Easy patient transport**

Stryker's Transport Stretcher features fifth wheel mobility for easier traction and cornering, and central-locking brakes for optimum stability.



Retractable fifth wheel provides improved mobility.



Broking system Metal braking system provides stability.



**Fold-down siderails** Fold-down siderails help ensure patient safety.

Option of **four-sided brake and steer controls** offer caregivers flexibility.

Optional **knee gatch** provides added patient comfort.

**Multiple bumper colour options** offers an opportunity to customise the stretcher for identification purposes.



Navy Grey Red

## stryker

# **Transport** Stretcher



### Standard features

- 20 cm (8 inch) Omni Surface casters
- Fold-down siderails on 66 cm (26 inch) litter (fold to foot)
- Fold-down siderails on 76 cm (30 inch) litter (fold to head)
- Four IV receptacles
- Central-locking brake system with dual-end activators
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Roller bumpers
- Stationary foot end litter
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side-mounted foot control hydraulics
- Hands-free Trendelenburg
- Wheel covers

### **Optional features and surface options**

- IsoFlex SE 76 cm (30") Surface width
- ComfortGel SE 66 or 76 cm (26"or 30") Surface width
- Ultra Comfort SE 76 cm (30") Surface width with choice of 8 or 10 cm (3" or 4") thickness, or 66 x 10 cm (26" x 4")
- Defibrillator tray

- Defibrillator tray/footboard/extender
- Footboard/chartholder
- Four-sided brake and steer controls
- Pop-up push handles
- IV poles

Warranty

- Permanent folding two- or three-stage at foot end
- Removable

- Patient restraints
- Serving/instrument tray
- Siderail pads
- Upright oxygen bottle holder
- Extended warranty available
- Crank knee gatch (model 748 only)

#### Medical

#### This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



01252005EN Rev J.3 DLS 03/2018 2017-16227 Copyright © 2018 Stryker stryker.com

Specifications	26-inch litter	30-inch litter
Model number	747	748
Overall length	211 cm (83")	211 cm (83")
Overall width	76 cm (30")	86 cm (34")
Weight capacity	227 kg (500 lb.)	227 kg (500 lb.)
Height range (to litter top)*		
High	91 cm (36")	91 cm (36")
Low	55 cm (21.5")	55 cm (21.5")
Litter positioning		
Backrest	0°- 90°	0°- 90°
Knee gatch		0°- 30°
Trend./Reverse Trend.	±18°	±18°
Patient surface	66 x 192 cm (26" x 75.5")	76 x 192 cm (30" x 75.5")
Siderails	33 x 140 cm (13" x 55")	33 x 140 cm (13" x 55")
Caster diameter	20 cm (8")	20 cm (8")

One year parts, labor and travel or two years parts. Warranty may vary by country. Contact your local Stryker Medical representative for additional information.