

Tru-Fit

Overbed tables

A true fit for every patient room

The sleek design of the award-winning Tru-Fit overbed table easily slides under beds, stretchers and mobile recliners.

Available in several laminate and base colours, our single-top and split-top overbed tables coordinate with any room décor.



U-shaped base

has a low profile to provide maximum versatility.

Infinite adjustments

between low and high heights.

High-impact

polystyrene cover protects the base and conceals the castors.

Dura-Edge

protects the table and contains spills.

Customisable

with vanity compartment and acrylic mirror.

Storage compartment


option provides additional storage space.



Single-top with vanity



Single-top without vanity



Split-top overbed table



Split-top with storage compartment



Tru-Fit

Overbed tables

Specifications

Model number	Single-top without vanity 3150-000-001	Single-top with vanity 3150-000-002	Split-top with vanity 3150-000-003
Height			
High	111 cm (43.75")	111 cm (43.75")	116 cm (45.5")
Low	69 cm (27")	69 cm (27")	74 cm (29")
Weight capacity*	34 kg (75 lb.)	34 kg (75 lb.)	34 kg (75 lb.)
Weight	20 kg (43 lb.)	20 kg (43 lb.)	28 kg (61 lb.)
Overall table top			
Length	79 cm (31")	79 cm (31")	79 cm (31")
Width	44 cm (17.5")	44 cm (17.5")	44 cm (17.5")
Base			
Height	5.6 cm (2.2")	5.6 cm (2.2")	5.6 cm (2.2")
Length	81 cm (32")	81 cm (32")	81 cm (32")
Width	54 cm (21.25")	54 cm (21.25")	54 cm (21.25")
Castor diameter	3.7 cm (1.5")	3.7 cm (1.5")	3.7 cm (1.5")
Storage compartment Part # 3151-050-001			
Weight capacity	18 kg (40 lb.)	18 kg (40 lb.)	19 kg (41 lb.)
Total unit with storage compartment			
Weight	31 kg (68 lb.)	43 kg (75 lb.)	39 kg (86 lb.)
Length/width	109 x 49 cm (42.75" x 17.5")	109 x 49 cm (42.75" x 17.5")	109 x 49 cm (42.75" x 17.5")

*11 kg safe working load

Key features

- 3.7 cm (1.5") diameter twin swivel castors
- Durable high-impact polystyrene base cover
- Height adjustment to 111 cm (43.75") (single-top only)
- Height adjustment to 116 cm (45.5") (split-top only)
- Infinite adjustments between low and high heights
- Anodised aluminium column
- Tabletop with Dura-Edge

Optional features

- Choose from six standard laminates and four standard base/Dura-Edge colours
- Storage compartment
- Secondary laminate top (split-top only)
- Vanity compartment with acrylic mirror

Warranty

- One year

For more information please contact your local sales representative.

Stryker Australia Pty Ltd
8 Herbert Street
St Leonards NSW 2065
61 2 9467 1000
www.stryker.com.au

Stryker New Zealand Limited
PO Box 17-136
Greenlane, Auckland, New Zealand
64 9 573 1890
www.stryker.com

Stryker reserves the right to change specifications without notice. Furniture is manufactured exclusively for Stryker Medical.

The information here is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.