



# Hands-free, wearable communication **featuring the Sync Badge**



## Hands-on patient care needs **hands-free communication**

Hands-on patient care happens anywhere the work is unpredictable, change is constant and teamwork is essential. Whether you need to call a code, call for equipment or call for help, you should always feel connected to the people you rely on. When your hands are on the patient – especially if you're wearing PPE – unlocking a smartphone, logging into an app and searching for a name in a directory might not be an option.

## Communication should **fit the way you work**

Inspired by nurses and designed for fast, reliable collaboration, the Sync Badge is a hands-free, wearable communication device that helps enhance workflow efficiency in dynamic, hands-on patient care environments.



### **Communicate as easily as you speak**

Say, "Okay Vocera" to wake up the Sync Badge and use voice commands to reach contacts by name, role, group or extension. For additional privacy, hold the Sync Badge to your ear, and it automatically switches to handset mode.

### **See and do the most important things first**

The Sync Badge's inbox organizes new messages, seen messages and calls, helping you to act promptly on communications and easily access previously viewed items.

### **Let the alarms and notifications you need come to you**

Receive alarms, critical lab results and STAT order notifications directly.\* Through integration with Stryker's Vocera Engage middleware, the Sync Badge can receive information from your EHR, nurse call, medical devices, scheduling systems and more.

### **Summon help at the press of a button**

Press the dedicated panic button in a moment of distress to quickly open a communication channel with your facility's response or security team. In a less emergent situation, you can use voice commands to broadcast to nearby colleagues for assistance.

### **Reach the entire care team seamlessly**

Communicate with team members across units and roles whether they use a hands-free, wearable device, a smartphone, a VoIP or landline phone, or a pager.

# Choose the right device for your role and workflow



## Sync Badge

Navigate intuitively and see what's most important with color coding and a three-section inbox.



## Smartbadge

Type text messages and view patient, event and care team context with notifications on the ample touchscreen.



## Minibadge

Triage urgency of events, call for help and more without needing to use your hands.

	Sync Badge	Smartbadge	Minibadge
<b>Dimensions</b>	98 mm x 38 mm x 17.5 mm	100 mm x 52 mm x 17 mm	62 mm x 52 mm x 17 mm
<b>Weight</b>	75 g	89 g	62 g
<b>Screen size</b>	Display with 1.3 in. color LCD 240 x 240 pixels	2.4" touchscreen with haptic touch call button	Not applicable
<b>Context displayed with notifications</b>	Patient and event context	Patient, event and care team context	Not applicable
<b>Ability to send text messages</b>	Not applicable	✓	Not applicable
<b>LED indicators for call status and presence visualization</b>	✓	✓	✓
<b>Dedicated panic button</b>	✓	✓	✓
<b>"Do not disturb" button</b>	✓	✓	✓
<b>Microphone array</b>	✓	✓	✓
<b>Hands-free speaker</b>	✓	✓	✓
<b>Dedicated handset mode receiver</b>	✓	✓	Use with headset or smartphone for handset mode
<b>USB-C charging port</b>	✓	✓	✓
<b>802.11ac dual-band wi-fi</b>	✓	✓	✓
<b>Supports shared devices</b>	✓	✓	✓
<b>Removable battery</b>	3.5 hours talk time, 45 hours standby	4.6 hours talk time, 37 hours standby	3.2 hours talk time, 41 hours standby
<b>8-bay charging stations</b>	✓	✓	✓
<b>Wearable accessories</b>	✓	✓	✓
<b>Integrated hardware, software and firmware</b>	✓	✓	✓
<b>System requirements: Consult with your Vocera account manager for specific feature compatibility.</b>	Consult with your Vocera account manager for specific feature compatibility	Vocera Voice Server 5.3.1 or later, Vocera Messaging Server 5.3.1 or later, Vocera Platform 6.1 or later	Vocera Voice Server 5.6 or later, Vocera Platform 6.5 or later

The names and characters presented throughout are entirely fictitious, any similarity to names or incidents is entirely coincidental.

\* US: Delivery of secondary alarms requires Engage Medical Device Alarm Notification (EMDAN), FDA 510(k)-cleared middleware, part of Stryker's Vocera Engage middleware. Capabilities are based on your facility's protocols and technology configuration. Canada: Delivery of secondary alarms requires Engage Medical Device Alarm Notification (EMDAN), part of Stryker's Vocera Engage middleware. Capabilities are based on your facility's protocols and technology configuration.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, Vocera Smartbadge, Vocera Minibadge, Vocera. All other trademarks are trademarks of their respective owners or holders. The absence of a product or service name or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.