

ST1-X Stretcher



Basic patient transport.

Stryker's ST1-X is your basic transport stretcher. Designed to satisfy your fundamental department needs, this stretcher includes fifth wheel mobility for easier steering and dual deck HPL x-ray for imaging needs.

Standard features

- Fifth wheel mobility
- Glideaway siderails
- Dual x-ray deck (HPL, full length)
- 4 integrated IV receptacles
- Central lock mechanism

Optional features

- Basic mattress 8 cm x 61 cm x 193 cm (3" x 24" x 75.5")
- Pop-up push handle options
- Restraint straps
- Hood base sticker (red, navy blue, gray)
- Upright oxygen bottle holder
- Defibrillator tray with chartholder
- Paper roll holder

Specifications

Model number	
ST1-X	6300-000-200
Overall length	217 cm (85.4")
Height	61 cm (24")
Overall width	78 cm (30.7")
Siderails up	78 cm (30.7")
Siderails down	73.5 cm (28.9")
Weight capacity	215kg (474 lbs.)
Trend./Reverse Trend.	±16°
Patient surface	62 ± 0.5 cm x 193 ± 0.5 cm (24.4" x 76")

Warranty

Warranty may vary by country. Contact your local Stryker Medical representative for additional information.

Medical

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker and ST1-X. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



PHMISC-BRO-1_19986EN Rev 1
SDL 09/2019
2019-19986

Copyright © 2019 Stryker
stryker.com