

stryker

AXS

Axial
Stability
Screw

**The Next Generation Screw
for Neuro Fixation**





AXS

(Axial Stability Screw)

Our latest iteration of the Universal Neuro Screw delivers on Stryker's commitment to advancing Neurotechnology. A brand new head design provides surgeons with the ease-of-use they desire with the performance they demand.

- Enhanced pick-up designed to help improve complications in loaded screwdriver hand-off from tech to surgeon, with optimized self-retention for reliable transport into OR field
- Self-centering feature facilitates off-axis insertion (fig.1)
- Low profile heights:
 - .45mm (3mm)
 - .55mm (4mm)
 - .65mm (5mm)



(fig.1)

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Stryker Australia
8 Herbert Street
St Leonards
NSW 2065 Australia
Ph: +61 2 9467 1000
Fax: +61 2 9467 1010

Stryker New Zealand
515 Mt. Wellington Highway
Auckland 1060 New Zealand
Ph: +64 9 573 1890
Fax: +64 9 573 1891

www.stryker.com
CMF-FL-2-SSP Rev. None
Copyright © 2016 Stryker
Printed in USA

CAT#	DESCRIPTION	(mm)
62-15035	Screwdriver blade long	
62-15036	Screwdriver blade short	
56-15034	Self-tapping screw	1.5x4mm
56-15035	Self-tapping screw	1.5x5mm
56-15036	Self-tapping screw	1.5x6mm
56-15933	Self-drilling screw	1.5x3mm
56-15934	Self-drilling screw	1.5x4mm
56-15935	Self-drilling screw	1.5x5mm
56-17334	Emergency screw	1.7x4mm