



2D and 3D **Universal Mesh**





Product portfolio

A legacy of performance

Our passion for quality, innovative technology, and delivering exceptional service spans more than half of a century. Crafted with a focus on fundamentals through relentless attention to detail, our products embody the spirit of Leibinger.



29-91002
Storage tray for
mesh instruments



01-01036
Mesh cutter



01-01038
Mesh clipper



2D Mesh

- **Countersink** - designed to reduce palpability
- Smooth edge can aid insertion
- MR conditional allows imaging after implantation
- **Variety of mesh types** - strength vs malleability
- **0.3mm Hybrid mesh** - 3x stronger than 0.3mm Dynamic mesh¹

3D Gold mesh

- **Countersink** - designed to reduce palpability
- Has stiffness equivalent to gold 2D mesh³
- Used for open defects up to 70mm and muscle attachment⁴
- Anatomically pre-formed mesh can minimize bending to fit patient
- MR conditional allows imaging after implantation

3D Silver mesh

- **Countersink** - designed to reduce palpability
- Increased stiffness for increased patient protection compared to 2D mesh³
- Used for open defects up to the size of the mesh and muscle attachment⁴
- Anatomically pre-formed mesh can minimize bending to fit patient
- MR conditional allows imaging after implantation

Mesh

Item number	Description	Size
1.2mm Dynamic mesh		
56-90322	Dynamic mesh - malleable - small	40x40x0.3mm
56-90324	Dynamic mesh - malleable - medium	90x90x0.3mm
56-90326	Dynamic mesh - malleable - large	120x120x0.3mm
56-90622	Dynamic mesh - standard - small	40x40x0.6mm
56-90624	Dynamic mesh - standard - medium	90x90x0.6mm
56-90626	Dynamic mesh - standard - large	120x120x0.6mm
56-90628	Dynamic mesh - standard - X-large	200x200x0.6mm
1.5/1.7mm Dynamic mesh²		
56-90312	Dynamic mesh - malleable - small	40x40x0.3mm
56-90314	Dynamic mesh - malleable - medium	90x90x0.3mm
56-90316	Dynamic mesh - malleable - large	120x120x0.3mm
56-90612	Dynamic mesh - standard - small	40x40x0.6mm
56-90614	Dynamic mesh - standard - medium	90x90x0.6mm
56-90616	Dynamic mesh - standard - large	120x120x0.6mm
56-90618	Dynamic mesh - standard - X-large	200x200x0.6mm
56-90814	Dynamic mesh - stiff - medium	90x90x0.8mm
56-90816	Dynamic mesh - stiff - large	120x120x0.8mm
56-90818	Dynamic mesh - stiff - X-large	200x200x0.8mm
1.2/1.5/1.7mm Hybrid mesh		
56-90342	Hybrid mesh - malleable - small	60x60x0.3mm
56-90344	Hybrid mesh - malleable - medium	90x90x0.3mm
56-90346	Hybrid mesh - malleable - large	120x120x0.3mm

Storage and instrumentation




2D mesh	
29-92001	2D mesh storage rack with lid
29-92002	2D mesh storage tray for implants with lid
29-92003	2D mesh silicone mat - small
29-92004	2D mesh silicone mat - large
29-91002	Storage tray for mesh instruments
01-01036	Mesh cutter coarse
01-01038	Mesh clipper

Dynamic mesh

Light blue

Gold

Green




Thickness: 0.3mm

Thickness: 0.6mm

Thickness: 0.8mm

Hybrid mesh



Thickness: 0.3mm

Item number	Description	Size
1.5/1.7mm Pre-form mesh		
56-90654	Pre-form mesh, standard - medium	90x90x0.6mm
56-90656	Pre-form mesh, standard - large	120x120x0.6mm
56-90658	Pre-form mesh, standard - XL	190x140x0.6mm

Item number	Description	Size
1.5/1.7mm Pre-form lateral mesh		
56-91064	Pre-form lateral mesh, stiff - left	175x120x80mm - 1.0mm
56-91062	Pre-form lateral mesh, stiff - right	175x120x80mm - 1.0mm

3D mesh	
29-93001	3D mesh storage rack with lid
29-93002	3D mesh storage tray for implants with lid
29-93003	3D mesh silicone mat
29-91002	Storage tray for mesh instruments
01-01036	Mesh cutter coarse
01-01038	Mesh clipper

Craniomaxillofacial

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets.

Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Distributed by:
Stryker Craniomaxillofacial
Kalamazoo, MI 49002 USA
t: 269 389 5346
toll free: 800 962 6558
f: 877 648 7114
stryker.com/cmf