

LIFEPAK® 15

monitor/defibrillator

Genuine accessories

LIFEPAK 15 monitor/defibrillator

Genuine accessories from **Stryker**

Ensure the safety of your staff and patients. Only Stryker accessories and disposables have been thoroughly tested with LIFEPAK products. Stick with the name you trust and the company that stands behind it. To order, contact your Stryker representative.



Power options



REDI-CHARGE® Base
AC power cord and adapter tray
not included.
11141-000115

**REDI-CHARGE LIFEPAK 15
Adapter Tray**
11140-000052



Mobile Battery Charger
Includes AC and DC power
cords, mounting bracket and
operating instructions.
11577-000011 (U.S.)
11577-000015 (United Kingdom)
11577-000012 (Europe)
11577-000017 (Australia)
11577-000014 (Switzerland)



**Rechargeable Lithium-ion
Battery**
5.7 ah, 10.8 volt.
21330-001176

AC Power Cords

11140-000015 (U.S.)
11140-000019 (Europe)
11140-000023 (South Africa)
11140-000046 (Italian)
11140-000017 (Australia)
11140-000021 (U.K.)
11140-000045 (Switzerland)



AC Power Adapter
Includes Right Angle
Cable (AC power cord
not included)
11140-000098

DC Power Adapter
Includes DC cable
and Right Angle Cable
11140-000074

AC Power Cord
11140-000015 (U.S.)
11140-000019 (Europe)
11140-000023 (South Africa)
11140-000046 (Italian)
11140-000017 (Australia)
11140-000021 (U.K.)
11140-000045 (Switzerland)

Power Attachment Kit
11577-000019



Extension Cable
For AC/DC Power Adapter.
11140-000080




**Replacement Right
Angle Power Cable**
For AC/DC Power Adapter.
11140-000081



**Replacement DC
Input Cable**
For DC Power Adapter.
11140-000084

ECG monitoring Accessories




12-Lead ECG Cable Trunk Cable with 4-Wire Limb Leads

11111-000018 (5ft, AHA)
11111-000019 (5ft, IEC)
11111-000020 (8ft, AHA)
11111-000021 (8ft, IEC)



12-Lead ECG Cable 6-Wire Precordial Attachment

11111-000022 (AHA)
11111-000023 (IEC)



4-wire Cable Comb
21300-008054 (10/pack)

6-wire Cable Comb
21330-008055 (10/pack)



3-Wire ECG Cable

11110-000029 (AHA)
11110-000030 (IEC)



5-Wire ECG Cable

4-wire limb plus 1 chest lead, labeled "V1" on the monitor reports.
11110-000066 (AHA)
11110-000067 (IEC)



Chart Recorder Paper

100mm x 22m
11240-000016 (2 rolls per box)



LIFE•PATCH® ECG Electrodes

Adult, pregelled
(Not available in the European Union)
11100-000001 (3/pack)
11100-000002 (4/pack)

Cleartrace® ECG Electrodes

Adult, clear tape
(Available in the European Union)
11100-000009 (3/pack)

Therapy delivery accessories

Hard paddles and electrode gel



Standard Hard Paddles (1 pair)

11130-000061 (English)
 11130-000062 (French)
 11130-000063 (German)
 11130-000064 (Spanish)
 11130-000065 (Italian)
 11130-000066 (Swedish)
 11130-000067 (Danish)
 11130-000068 (Portuguese)
 11130-000069 (Norwegian)
 11130-000070 (Dutch)
 11130-000071 (Finnish)
 11130-000072 (Polish)
 11130-000073 (Japanese)



Pediatric Paddle Adapter

11133-000007 (1 adapter)



SIGNAGEL® Electrode Gel

For use with hard paddles.
 Recommended for ECG, defibrillation,
 biofeedback and EMG.
 21300-005847 (8.5 oz)



Internal Paddles (Requires Internal Paddle Handles and Internal Paddles Adapter Cable)

1 in. size	2.5 in. size
11131-000010 (1 pair, 6.25 in. shaft)	11131-000013 (1 pair, 5.5 in. shaft)
	11131-000019 (1 pair, 8.5 in. shaft)
1.5 in. size	3.5 in. size
11131-000011 (1 pair, 6 in. shaft)	11131-000014 (1 pair, 5 in. shaft)
11131-000021 (1 pair, 9 in. shaft)	11131-000023 (1 pair, 8 in. shaft)
11131-000024 (1 pair, 14 in. shaft)	

2 in. size
 11131-000012 (1 pair, 5.75 in. shaft)
 11131-000022 (1 pair, 8.75 in. shaft)



Internal Paddle Handles with Discharge Control

(For use with the Internal Paddles Adapter Cable)
 11131-000001 (1 pair)



Internal Paddles Adapter Cable
 (For use with Internal Paddle Handles)
 11998-000326

Therapy delivery accessories

EDGE System™ Electrodes for pacing/defibrillation/ECG with QUIK-COMBO® Connector

18-month minimum shelf life remaining at time of shipment from Physio-Control except where noted.



EDGE System Electrodes with QUIK-COMBO Connector
24" leadwire length
11996-000091



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector
24" leadwire length
11996-000090



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK® Preconnect System
42" leadwire length
11996-000017



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector
For use only with manual monitor/defibrillators; 12 month minimum shelf life at time of shipment 24" leadwire length.
11996-000093



QUIK-COMBO Therapy Cable
With convenient TRUE-LOCK® Cable Connector.
11113-000004 (8 ft)

NIBP monitoring accessories

NIBP hoses



NIBP Tubing

21300-008159 (6ft)
21300-008147 (9ft)
21300-008146 (12ft)



NIBP Tubing, Coiled

21300-008148 (2-9ft)

NIBP cuffs



Reusable Cuff

X-Large Adult 35 - 44 cm 11160-000019	Pediatric 13 - 20 cm 11160-000013
Large Adult 32 - 42 cm 11160-000017	Infant 8 - 14 cm 11160-000011
Adult 2 6 - 35 cm 11160-000015	



Single Patient Use Cuff

X-Large Adult 35 - 44 cm 11160-000020	Pediatric 13 - 20 cm 11160-000014
Large Adult 32 - 42 cm 11160-000018	Infant 8 - 14 cm 11160-000012
Adult 26 - 35 cm 11160-000016	

Pulse oximetry monitoring accessories

Masimo SET® RC patient cables



RC Patient Cable

For use with M-LNCS and Rainbow Patient Sensors

11171-000069 (1ft)

11171-000037 (4ft)

11171-000038 (12ft)



EMS RC Patient Cable

11171-000082 (4 ft)

Masimo SET RC patient cable compatible SpO₂ sensors



M-Incs Reusable Sensor

11171-000046 (Ad)


11171-000047 (Ped)



M-LNCS Adhesive Sensors (20/box)

11171-000039 (Ad)

11171-000040 (Ped)




M-LNCS Adhesive Sensors (20/box)

11171-000043 (Neo/Pt)

11171-000042 (Neo/Ad)

11171-000041 (Inf)

Masimo SET RC patient cable compatible Rainbow® SpO₂, SpCO, SpMet sensors



Rainbow Reusable Sensor

11171-000049 (Ad)

11171-000050 (Ped)



Rainbow Adhesive Sensor (10/box)

11996-000339 (Ad)

11996-000340 (Ped)



Rainbow Adhesive Sensor (10/box)

11996-000342 (Inf)

11996-000341 (Neo/Ad)

Masimo SET LNC patient cables



Red LNC Patient Cable

For Use with LNCS Patient Sensors

11996-000323 (4ft)

11996-000324 (10ft)

11996-000325 (14ft)

Masimo SET LNC patient cable compatible SpO₂ sensors



LNCS® Reusable Sensor


11171-000017 (Ad)

11171-000018 (Ped)



LNCS Reusable Soft Sensor

11171-000052 (Ad)



LNCS Adhesive Sensor (20/box)

11171-000019 (Ad)

11171-000020 (Ped)



LNCS Adhesive Sensor (20/box)

11171-000029 (Neo/Pt)

11171-000028 (Neo/Ad)

11171-000031 (Inf)

Pulse oximetry monitoring accessories

Direct connect SpO₂ only patient sensors



Adult Reusable Direct Connect Sensor

11996-000331 (3ft)
11996-000332 (12ft)



Pediatric Reusable Direct Connect Sensor

11996-000333 (3ft)
11996-000334 (12ft)



Adult Reusable Soft Direct Connect Sensor

11171-000053 (8ft)

Direct connect Rainbow SpO₂, SpCO, SpMet patient sensors



Adult Rainbow Direct Connect Reusable Sensor


11996-000335 (3ft)
11171-000032 (8ft)
11996-000336 (12ft)



Pediatric Rainbow Direct Connect Reusable Sensor

11996-000337 (3ft)
11171-000033 (8ft)
11996-000338 (12ft)

Additional Masimo accessories



Reusable Ambient Light Shield

11171-000054 (5/bag)



Disposable Ambient Light Shield

11171-000055 (10/bag)

Pulse oximetry monitoring accessories

Masimo to Nellcor™ adapter



Red MNC Cable

Connects LIFEPAK 15 to Nellcor patient sensor

11996-000365 (4ft)

11996-000366 (10ft)

Masimo to Nellcor adapter compatible SpO₂ sensors



DURASENSOR Reusable Clip

11996-000060 (Ad)

DURA-Y Multisite Reusable Sensor

11996-000106 (> 1 kg)



Oxiband Reusable Sensor

Includes 50 disposable adhesive sensors.

11996-000061 (Ad/Neo)

11996-000062 (Ped/Inf)

Disposable Adhesive Bandage Wrap (100/pk)

Not Available in Canada.

11996-000048 (Ad/Neo)

11996-000049 (Ped/Inf)

Temperature monitoring

Temperature cables



Temperature Adapter Cable

11140-000079 (10ft)

11140-000078 (5ft)

Temperature sensors



Esophageal-Rectal Sensor

11996-000360 (9FR, 20/box)



Skin Probe Sensor

11996-000359 (20/box)

End-Tidal CO₂ (EtCO₂) monitoring accessories

Filterlines for intubated patients

Single patient use



FilterLine® SET

Key Applications: OR, EMS, ED, Rapid Response Teams, Transport.

Adult/Pediatric

11996-000081 (25/pk, 200 cm)

11996-000164 (25/pk, 400 cm)



FilterLine H SET

Key Applications: Critical care, Humidified Environments.

Adult/Pediatric

11996-000080 (25/pk, 200 cm)

Infant/Neonatal

11996-000001 (25/pk, 200 cm)

Non-intubated filterlines

Single patient use



Smart CapnoLine® Plus

Key Applications: Procedural Sedation, Upper GI Procedures, MAC, RMS, BD, Rapid Response Teams.

Adult with O₂

11996-000163 (25/pk, 200 cm)

11996-000167 (100/pk, 200 cm)

11996-000165 (25/pk, 400 cm)

Adult without O₂

11996-000162 (25/pk, 200 cm)

11996-000166 (100/pk, 200 cm)



Smart CapnoLine

Key Applications: Procedural Sedation, SMS, SD, Rapid Response Teams.

Pediatric with O₂

11996-000128 (25/pk, 200 cm)

Pediatric without O₂

11996-000120 (25/pk, 200 cm)

LIFEPAK 15 monitor/defibrillator

Cases & mounting options



Standard Carrying Case

Includes right pouch and left pouch.
11577-000002

Shoulder Strap

11577-000001



Bed Connector

For use in hospital only.
11996-000374

Top Pouch

Storage for sensors and electrodes;
insert in place of standard paddles.
11220-000028

Back Pouch

Ideal for additional accessory storage.
11260-000039

Communication accessories



LIFEPAK Monitor to PC Cable

For connecting LIFEPAK 12 or
LIFEPAK 15 monitor/defibrillator
to PC.

11230-000020 (Serial)
11996-000369 (USB)



Titan III Wireless Gateway

For transmitting data from
LIFEPAK 15 to the LIFENET®
System or CODE-STAT™ data
review software. Requires
existing wireless network.

21996-000115 (Wireless Gateway)
21996-000114 (Wireless, 3G + Voice
Gateway)
21996-000116 (Wireless 3G Gateway)



3G Gateway

For transmitting data from
LIFEPAK 15 to the LIFENET
System or CODE-STAT data review
software. Requires data plan.

21996-000087

Testers & training materials




Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias, and bradycardias, as well as ST segment and T wave abnormalities.

11996-000311 (12-Lead)

11996-000310 (3-Lead)



Defibrillator Checker

Tests integrity of energy delivery through standard hard paddles. Neon light indicates energy has been delivered.

11998-000060



Test Load

Use to perform therapy cable performance checks. Connects to the QUIK-COMBO therapy cable on the defibrillator.

21330-001368 (Brazilian, Hungarian, Portuguese, Polish, Romanian)

21330-001440 (Chinese, Czech, Japanese, Korean, Russian)

21330-001366 (Danish, Finnish, Norwegian, Swedish)

21330-001367 (Dutch, French Canadian, German, Intl English, Italian, Spanish)

21330-001365 (English)

Training tools



Inservice Video:

21330-001357 (English - NTSC DVD)

Operating Instructions:

26500-003978 (English)
26500-003979 (Intl English)
26500-003980 (German)
26500-003981 (Italian)
26500-003982 (French)
26500-003983 (Dutch)
26500-003984 (Spanish)
26500-003985 (Iberian Portuguese)
26500-003986 (Brazilian Portuguese)
26500-003987 (Swedish)
26500-003988 (Danish)
26500-003989 (Finnish)
26500-003990 (Norwegian)
26500-003991 (Polish)
26500-003992 (Hungarian)
26500-003993 (Czech)
26500-003994 (Russian)
26500-003995 (Korean)
26500-003996 (Japanese)

**For further information, please contact your
Stryker representative or visit our website at
strykeremergencycare.com**

Emergency Care

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CODE-STAT, EDGE System, LIFEPAK, REDI-PAK, QUIK-COMBO, REDI-CHARGE, Stryker, TRUE-LOCK, Masimo, the Radical logo, Rainbow and SET are registered trademarks of Masimo Corporation. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

GDR 3301060_L
Copyright © 2019 Stryker

Distributed by:

Stryker
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
strykeremergencycare.com

Stryker European
Operations B.V.
Herikerbergweg 110
1101 CM Amsterdam
Netherlands
Tel +31 (0)43 3620008
Fax +31 (0)43 3632001

Distributed by:

Stryker UK Ltd
Stryker House
Hambridge Road
Newbury, Berkshire
RG14 5AW
United Kingdom

Distributed by:

Stryker Australia Pty Ltd
8 Herbert Street
St Leonards NSW 2065
Australia
Toll Free Tel 1800 987 982
Toll Free Fax 1800 890 892

Stryker New Zealand Limited
515 Mt Wellington Highway
Mt Wellington, Auckland 1060
New Zealand
Tel +61 2 9415 5900
Fax +61 2 80767672