

CODE-STAT™ 11 Data Review Software

System Recommendations

The hardware and software described below include the minimal and recommended system capabilities needed to set up and run CODE-STAT 11 data review software.



Processor	Minimal: Pentium® IV-class 2.4 GHz Recommended: Pentium IV-class 3.0 GHz
Operating System	Windows® 7 with Service Pack 1 (32-bit and 64-bit) Windows 10 (64-bit) Windows Server 2008 R2 with Service Pack 1 (64-bit) Windows Server 2012 R2 (64-bit) Latest Microsoft® updates recommended
Database	Microsoft SQL Server 2014 (included) Customer-supplied SQL Server 2008 R2 SP2 or later (must have existing licenses and permissions)
Software	Microsoft .NET Framework 4.6.2 Client Profile required Microsoft .NET Framework 3.5 is required to install CODE-STAT with a local database Windows Media Player is required to play standalone audio and video reports
Memory	Minimal for 32-bit operating systems: 1 GB RAM Minimal for 64-bit operating systems: 2 GB RAM Recommended: 4 GB RAM
Hard Disk	A minimum of 2 GB free space
Sound	Windows-compatible sound card and speakers (required to play back audio reports that are associated with continuous ECG)
Screen Resolution	Minimal: 1024 x 768 pixels

Physio-Control is now part of Stryker.

©2018 Physio-Control, Inc. Redmond, WA. All names herein are trademarks or registered trademarks of their respective owners.
GDR 3336982_A