


Better view. Better care.

See clearly with the next-generation McGRATH™ MAC video laryngoscope.

You need a simple and convenient intubation solution. The enhanced McGRATH MAC video laryngoscope helps your team:

- Increase visibility with improved optics and a brighter, warmer LED light.
- Gain peace of mind during any procedure with a more durable handle.
- Automatically save battery life during idle time with auto-off feature.



2x light spread
3x brighter
4x resolution¹



	Legacy design	Enhanced design
Visualization		
Light color	Bright, white light	Natural, warm light
Light brightness	Standard light source	3x brighter light source ²
White balance	Fixed	Automatically adjusts
Camera detail	Standard resolution	Enhanced resolution
Durability		
McGRATH MAC handle	Standard screen hinge	Reinforced screen hinge
IP rating	IPX7	IPX8
Battery management		
	Manual power on/off feature	Auto-off feature
Ordering information		
SKU number	11996-000393	301-000-000



Contact your local Stryker account manager for details.

1. Compared to the previous version of the McGRATH MAC video laryngoscope.
 2. Based on internal test report #RE001500188, V&V plan.

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate Stryker's product offerings. A healthcare professional must always refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

All names herein are trademarks or registered trademarks of their respective owners. McGRATH™ MAC is a registered trademark of Medtronic. Not all products and services are available in all countries.

Distributed by:
 Stryker
 11811 Willows Road NE
 Redmond, WA, 98052 U.S.A.
 Toll free 800 442 1142
 strykeremergencycare.com

Distributed in Canada by:
 Stryker Canada
 2 Medicorum Place
 Waterdown, Ontario
 L8B 1W2
 Canada
 Toll free 800 668 8323