


# Performance-LOAD<sup>®</sup>

## Manual Cot Fastener



3800 E. Centre Avenue  
 Portage, MI 49002 USA  
 t: 269 329 2100  
 toll free: 800 327 0770  
[ems.stryker.com](http://ems.stryker.com)

### Key Features

- **SAE J3027 compliant** to maximize occupant safety in the event of an ambulance accident
- **Guided loading and unloading** help minimize cot drifting during loading and unloading to improve operator efficiencies
- **Adaptable design** with the ability for future modification to support Stryker Power-LOAD fastener system
- **Mass Casualty compliant** with the inclusion of wall or floor mount fasteners and appropriate hardware
- **Power washable and easily removable** for maximum decontamination and lowered potential for biohazards

### Optional Features

- **Inductive charging** provides optimal SMRT battery management and charging when fastened (no cables or connectors required)

### Warranty

- Two-year parts and labor
- Lifetime on all welds\*
- \*7-year service life

**Extended warranties available.**

### Specifications\*

<b>Model Number</b>	6392
<b>Length</b>	
Overall Length	70 in (178 cm)
Length with Mass Casualty inclusion	80 in (203 cm)
<b>Width</b>	
Foot End	9 in (23 cm)
Head End	19 in (48 cm)
<b>Weight</b>	
Total Weight	63 lb (28.5 kg)
Floor Plate Assembly	17 lb (7.5 kg)
Transfer Assembly	46 lb (21 kg)
<b>Height</b>	6 in (15 cm)
<b>Minimum Operator Required</b>	
Occupied Cot	2
Unoccupied Cot	1
<b>Recommended Loading Height</b>	22 in to 36 in (56 cm to 91 cm)

Meets dynamic crash standards with Power-PRO XT (AS/NZS-4535, BS EN-1789 with X-restraints and SAE J3027) with X-restraints and Performance-PRO XT (SAE J3027 with X-restraints).

Performance-LOAD is rated IP26 with inductive charging option

Stryker reserves the right to change specifications without notice.

Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

The Performance-LOAD cot fastener system is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Performance-LOAD, Power-LOAD, Stryker**, . All other trademarks are trademarks of their respective owners or holder.