

# Power-LOAD<sup>®</sup>

## Cot Fastener System



3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770  
ems.stryker.com

### Standard Features

- Lifting arms
- Head end lock LED indicators
- Control panel
- Cot release handles
- Linear transfer system
- Inductive charging
- Operation guide

### Optional Features

- Wheel guide
- Mass casualty floor mount assembly
- Mass casualty wall mount assembly

### Specifications

<b>Model Number</b>	<b>6390</b>
<b>Length</b>	
Overall Length	95 in (241 cm)
Minimum Length	89.5 in (228 cm)
<b>Width</b>	24.5 in (62 cm)
<b>Weight</b>	
Total Weight	211.5 lb (96.5 kg)
Floor Plate Assembly	16.5 lb (7.5 kg)
Anchor Assembly	23 lb (10.5 kg)
Transfer Assembly	67 lb (30.5 kg)
Trolley Assembly	105 lb (48 kg)
<b>Maximum Lift Capacity*</b>	700 lb (318 kg)
<b>Minimum Operator Required</b>	
Occupied Cot	2
Unoccupied Cot	1
<b>Recommended Loading Height</b>	22 in to 36 in (56 cm to 91 cm)
<b>Battery</b>	12 VDC, 5 Ah Lead Acid Battery (6390-001-468)

### Warranty

- One-year parts, labor, and travel

### Extended warranties available.

7-year expected service life.

Maintenance agreements are available for service beyond its expected service life.

\*Maximum weight capacity represents patient weight and accessories. Safe working load of 870 lb (395 kg) represents the sum of the cot total weight and patient.

Stryker reserves the right to change specifications without notice.

The Power-LOAD cot fastener system is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822.

The yellow and black color scheme is a registered trademark of Stryker Corporation.

Meets dynamic crash standards with Power-PRO™ XT cot fastening system (AS/NZS-4535, BS EN-1789 and SAE J3027 with X-restraints)

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Power-LOAD, Stryker**, . All other trademarks are trademarks of their respective owners or holder.

#### Certifications

IPX6 IEC-60601-1 BS EN-1789  
AS/NZS-4535 KKK-A-1822

Printed in U.S.A.

Copyright © 2018 Stryker

Mkt Lit-539 08 JUN 2011 Rev E