

Mini-PRO

Pediatric Restraint



Secure. Seamless. Standardized.

When it comes to caring for your patients, compromising isn't an option. Mini-PRO, Stryker's first pediatric restraint, allows you to deliver the necessary care to your patient while providing secure transport to the hospital.



Designed for transportation safety

- **Dynamically crash tested** to the pulses of SAE J3027, EN 1789 and AS/NZS 4535.
- **Color-coded restraints and D-rings** are designed to allow care providers to assemble Mini-PRO on the cot quickly.
- **Provide calm** in the chaos with co-transport capabilities for patients from 4 - 22 lbs (1.8 kg - 9.98 kg).
- **Deliver the necessary care** without compromising security - patient chest clasp can be removed to gain access to the patient's chest while remaining crash tested.
- **One size fits all solution:** 4 - 108lbs weight accommodation (1.8 kg - 49 kg).
- **Lateral head support**
- **Compatible** with Stryker's Power-PRO 2, Power-PRO XT, Performance-PRO and Power-PRO TL with the D-ring X-restraints.
- **Cleaning and disinfection:** Mini-PRO can be cleaned and disinfected following the same methods and chemicals listed in the Power-PRO manuals for restraints and mattress.*
- **7 year** expected service life

Equip your fleet with an intuitive pediatric restraint while keeping everyone in the back of the ambulance secure during transport.



Carrying case dimensions - 17.5 in X 12 in X 7 in
(44.45 cm x 30.48 cm x 17.78 cm)

Quick-start guide included

Mini-PRO is not yet, and may not be, offered for sale or lease, or sold or leased, until commercially available.

*Excludes power washing

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a patient and must refer to the instructions for use before using any Stryker product. Stryker Corporation or its divisions or other corporate affiliated entities own, use, or have applied for the following trademarks or service marks: Mini-PRO, Performance-PRO, Power-PRO and Stryker. All other trademarks are trademarks of their respective owners or holders.

04/2026
EG-PP2-FACT-3093301_REV-0_en_us | Copyright © 2026 Stryker



3800 E. Centre Avenue
Portage, MI 49002 U.S.A.
t: 269 329 2100
toll free: 800 327 0770

stryker.com