

stryker

Stair-PRO

stair chair



Stair-PRO

stair chair

Delivering proper patient care poses unique challenges. Patient handling in stairways is a challenge you may face on any given call. The strain on your back from lifting and lowering patients down stairs can add up over time. With 1 in 4 EMS workers suffering a career-ending back injury within the first four years of employment in the US alone*, we are striving to offer options that help reduce this risk.

Each patient handled in a stairway can pose a risk for caregiver injury. Stair-PRO puts the power to help reduce that risk back in your hands. Stryker's Stair-PRO stair chair – **it's a powerful thing.**



Compact folded depth

Easy-to-fold functionality allows the chair to fold to a compact size for convenient storage in existing vehicle compartments.

Optional accessories:



O₂ Holder



Foot Support



IV Clip

Key features:



Specifications¹

• Models	6252	
• Height	95.25 cm	37.5 in
• Width	52 cm	20.5 in
• Depth²	71 cm	28 in
• Folded Depth	20 cm	8 in
• Weight	14.9 kg	32.8 lb
• Maximum Load³	228 kg	500 lb

¹ Dimensions are measured from the outermost edges of the main frame. Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

² Depth dimensions are measured with extendable handles retracted.

³ Maximum load capacity is total weight distributed in accordance to basic human anatomy. EMTs must consider the weight of the patient, equipment and accessories when determining the total weight on the product.

⁴ Patient restraint set includes two cross-chest restraints and one ankle restraint. Cross-chest restraints can also be used in a chest and lap restraint configuration.

In-service video included with every order.

Stryker reserves the right to change specifications without notice.

References

* Sanders, Mick J. (2011) Mosby's Paramedic Textbook (4th ed., p. 36)

Technical Support

Stryker Technical Support comprises a team of professionals available to help with your Stair-PRO needs.

ProCare Service from Stryker

Every day, you count on your medical equipment to perform at its best. With ProCare Services, our highly trained technicians help ensure your equipment is ready to perform when needed. Our goal is to make life easier for your staff and clinicians with service plans designed to help you get the most out of your investment.

Flex Financial Program

Our financial programs provide a range of smart alternatives designed to fit your organization's needs. We offer flexibility beyond a cash purchase with payment structures customizable to meet budgetary needs and help to build long-term financial stability. Specialised finance solutions for municipalities and other government agency businesses are also available.

Program availability will vary depending on country and product portfolio, please check with your local sales representative for full details and availability.

Medical EMS

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.

Stryker's Regional Headquarters
Amsterdam Herikerbergweg 110
1101 CM, Amsterdam
The Netherlands



Mkt Lit-30501EN Rev D
DLS 01/2018
2017-15526

Copyright © 2018 Stryker
stryker.com