

MX-PRO[®]

Bariatric Transport Cot



3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770
www.ems.stryker.com

Standard Features

- Transport at load height capability
- Positive action height adjustment
- Easy-to-use release handle design
- Color-coded controls
- Lap belt extensions
- High visibility powder-coated frame
- Lightweight, aluminum construction
- Scientifically optimized lift bar design
- Lower lifting bar
- Seven height positions
- Integrated bumper system
- Lift-capable safety bar
- Perforated litter surface
- One-hand release breakaway head section with safety bar
- Floor-mounted safety hook
- One-hand release, fold-down siderails
- One-hand release, infinite positioning, pneumatically assisted backrest
- Oversized wheels with sealed caster and wheel bearings
- Reflective labeling
- Sealed flat mattress
- Shock positioning
- Two lap belts and one four-point shoulder restraint

Optional Features

- Heavy duty two or three-stage IV poles (patient right or left)
- Base storage tray
- Nonpocketed head end storage
- Height limit kit¹
- Permanent or removable oxygen bottle holders (foot end)
- Rigid head and foot end push/pull handles
- Single or dual wheel lock
- X-frame guards
- Side lift handles
- Tow package

Specifications

Model Number	6083
Height²	
Position 1	13.5 in. (34.3 cm)
Position 2	24.5 in. (62.2 cm)
Position 3	27.5 in. (69.9 cm)
Position 4	30 in. (76.2 cm)
Position 5	32.5 in. (82.6 cm)
Position 6	34.5 in. (87.6 cm)
Position 7	37 in. (94 cm)
Length	
Standard	80.5 in. (204.5 cm)
Minimum	62 in. (157.5 cm)
Width	29 in. (73.7 cm)
Weight³	111 lb (50.3 kg)
Wheels	
Diameter	6 in. (15.2 cm)
Width	2 in. (5.1 cm)
Articulation	
Backrest	2– 73°
Shock Position	+14°
Maximum Weight Capacity⁴	850/1600 lb (385.6/725.7 kg)
Minimum Operator Required	
Occupied Cot	2
Unoccupied Cot	1
Recommended Fastener System	
Floor Mount	Model 6370, 6377 6378
Wall Mount	Model 6371
Recommended Loading Height	Up to 32 in (81.3 cm)

Warranty

- One year parts, labor or two years parts
- One year parts and labor or fifteen months parts and labor
- Two year parts only
- Lifetime on all welds

Extended warranties available.

5-year service life.

¹ Height limit kit recommended for load heights less than 28 inches (71 cm). High height kit recommended for load heights more than 32 inches (81 cm).

² Height measured from bottom of mattress at seat section to ground level.

³ Cot is weighed without mattress and restraints.

⁴ 1600 lb capacity at the lowest position.

Stryker reserves the right to change specifications without notice.

The MX-PRO Bariatric Transport is designed to be compatible with competitive cot fastener systems.

Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

The yellow and black color scheme is a proprietary trademark of the Stryker Corporation.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **MX-Pro, Stryker**, . All other trademarks are trademarks of their respective owners or holder.