

WEEE recycling passport

LIFEPAK® 15, 20e defibrillator / monitor and accessories

Product Information

Product Name	LIFEPAK® 15, 20e defibrillator/monitor accessories
--------------	--

Producer/Manufacturer

Name of the Company	Stryker Medical
Address	11811 Willows Road NE, Redmond, WA 98052 USA

Recycling information (device)

EEE components	Present in the device		Materials, location and important information
	Yes	No	
printed circuit boards	x		lead in PCBA components and some solders under RoHS exemptions
batteries	X		lithium ion
coin cell battery	X		lithium battery on main PCBA
electrolyte Capacitor	x		note: the energy stored in the electrolytic capacitor is lethal when charged.
liquid crystal display	x		front of device, surface is less than 100cm ²
plastics	x		PET, PC+ABS, nylon, PTFE polypropylene, liquid crystal polymer, polycarbonate film, silicone, polyimide, polyolefin, PC+PET, polyurethane, PVC thermoplastic, thermoplastic PBT, acetal copolymer, glass fiber reinforced polyphenylene sulfide
plastics containing brominated flame retardants		x	N/A
connectors, fasteners, standoffs	x		inside of device, lead in brass RoHS exemptions
mercury		x	N/A
toner cartridges		x	N/A
asbestos		x	N/A
radioactive substances		x	N/A
chlorofluorocarbons (CFC), hydrochlorofluorocarbons		x	N/A

Emergency Care

11811 Willows Road NE, Redmond, WA 98052 USA | P +1 425 867 4000 | Toll-free +1 800 442 1142 | stryker.com

(HCFC) or hydrofluorocarbons (HFC), hydrocarbons (HC),			
gas discharge lamps		x	N/A
cathode ray tubes		x	N/A
polychlorinated biphenyls (PCB) containing capacitors		x	N/A

Recycling information (accessories)

Recycling/Hazardous Materials	Present in the Device		Important Information
	Yes	No	
printed circuit boards	x		in connectors
batteries		x	N/A
electrolyte capacitor		x	N/A
external electric cables	x		
liquid crystal display		x	N/A
plastics	x		PET, PC+ABS, nylon, PTFE polypropylene, liquid crystal polymer, polycarbonate film, silicone, polyimide, polyolefin, PC+PET, polyurethane, PVC thermoplastic, thermoplastic PBT, acetal copolymer, glass fiber reinforced polyphenylene sulfide
plastics containing brominated flame retardants		x	N/A
connectors, fasteners, standoffs	x		lead solder, lead in brass RoHS exemptions
mercury		X	N/A
toner cartridges		X	N/A
asbestos		X	N/A
radioactive substances		X	N/A
chlorofluorocarbons (CFC), hydrochlorofluorocarbons (HCFC) or hydrofluorocarbons (HFC), hydrocarbons (HC),		X	N/A
gas discharge lamps		X	N/A
cathode ray tubes		X	N/A
polychlorinated biphenyls (PCB) containing capacitors		X	N/A

Emergency Care

11811 Willows Road NE, Redmond, WA 98052 USA | P +1 425 867 4000 | Toll-free +1 800 442 1142 | stryker.com