

**REACH SVHC declaration**  
**LIFEPAK® 35 defibrillator / monitor and accessories**

**Product Information**

|              |   |
|--------------|---|
| Product Name | LIFEPAK® 35 defibrillator/monitor accessories |
|--------------|---|

**Producer/Manufacturer**

|                     |  |
|---------------------|--|
| Name of the Company | Stryker Medical                              |
| Address             | 11811 Willows Road NE, Redmond, WA 98052 USA |

**EU REACH SVHC Communication**

The Registration, Evaluation, Authorization, and Restriction of Chemicals (REACH), a European Union regulation, requires Stryker Medical Redmond to provide chemical content information for Substances of Very High Concern (SVHC) if they are present above 0.1% of the article weight. The SVHC list is updated on a regular basis.

**EU REACH Article 33 SVHC**

| Product                                   | Component                      | Substance | Cas Number |
|---|--------------------------------|-----------|------------|
| LIFEPAK® 35                               | Printed circuit board assembly | Lead (Pb) | 7439-92-1  |
|   | Electrical contacts            | Lead (Pb) | 7439-92-1  |
|   | Stand-offs                     | Lead (Pb) | 7439-92-1  |
|   | Keypad                         | Lead (Pb) | 7439-92-1  |
| LIFEPAK® 35 accessories<br>ECG wire cable | Electrical contacts            | Lead (Pb) | 7439-92-1  |

**Emergency Care**

11811 Willows Road NE, Redmond, WA 98052 USA | P +1 425 867 4000 | Toll-free +1 800 442 1142 | [stryker.com](http://stryker.com)