



Document number:

M0000032174

Name:

FM: Digital Product Passport

Revision:

AA.0

Reference:

M0000020673

Form

Digital Product Passport

Medical

Product Family: PEC – Alvarium™ Battery

Xpedition (model 6257)

Power-PRO 2 (model 6507)

Battery Model Number(s): 650700080401





Document number: M0000032174
Name: FM: Digital Product Passport
Revision: AA.0
Reference: M0000020673

Form

Batteries and Waste Batteries Regulation Compliance

General Battery Information

Battery Manufacturer & address	Design: Stryker Medical 3800 E. Centre Ave Portage, MI 49007 Manufacturer: Carretera Fresnillo-Plateros KM2 SN Fresnillo, Zacatecas, Mexico, CP99056
Battery Category	Industrial with capacity <2kWh <i>Annex VI -2: the battery category and information in accordance with Article 38(6): Manufacturers shall ensure that batteries which they place on the market bear a model identification and batch or serial number, or product number or another element allowing their identification. Where the size or nature of the battery does not allow it, the required information shall be provided on the packaging or in a document accompanying the battery.</i>
Place of Manufacturer	Carretera Fresnillo-Plateros KM2 SN Fresnillo, Zacatecas, Mexico, CP99056
Date of Manufacture	See product label
Weight	3.8 lbs. (1.7 kg)
Capacity	25.6V DC
Chemistry	Lithium-iron phosphate
Hazardous substances (other than mercury, cadmium, lead)	Lead, lead compounds, diborane trioxide, bisphenol A (BPA), antimony oxide (antimony trioxide)
Usable Extinguishing Agent	Dry chemical-type extinguishers or CO2 extinguishers
Critical Raw Materials Present in Battery in a Concentration Greater than 0.1% w/w.	Lithium iron phosphate (LiFePO4)- 25-30% Iron (FE)-22-30%
Instructions and Safety Information	See product manual; reference part numbers: 6507-009-001 and 6257-109-001
Storage and Transport Conditions	-38 degrees F (-34 degrees C) temperature 10%-93% relative humidity 500-1060 h Pa atmospheric pressure