



Document number:

M0000032174

Name:

FM: Digital Product Passport

Revision:

AA.0

Reference:

M0000020673

Form

Digital Product Passport

Medical

Product Family: Sealed Lead Battery

Power-LOAD (model 6390)

Battery Model Number(s): 6390-001-486





Document number: M0000032174
Name: FM: Digital Product Passport
Revision: AA.0
Reference: M0000020673

Form

Batteries and Waste Batteries Regulation Compliance

General Battery Information

Battery Manufacturer & address	EnerSys EnergyProductsInc. 617 N. Ridgeview Drive Warrensburg, MO 64093-9301
Battery Category	Industrial with capacity <2kWh <i>Annex VI -2: the battery category and information in accordance with Article 38(6): Manufacturers shall ensure that batteries which they place on the market bear a model identification and batch or serial number, or product number or another element allowing their identification. Where the size or nature of the battery does not allow it, the required information shall be provided on the packaging or in a document accompanying the battery.</i>
Place of Manufacturer	617 N. Ridgeview Drive Warrensburg, MO 64093-9301 Supplier: Midwest Safety Products 4929 East Paris Avenue, SE Grand Rapids, MI 49512
Date of Manufacture	See Product Label
Weight	5.9 lbs/ 2.7 kg
Capacity	12V, 5 AH LEAD Acid Battery
Chemistry	Sealed Lead Acid Battery, VRLA Battery
Hazardous substances (other than mercury, cadmium, lead)	Sulfuric Acid, Lead Compounds
Usable Extinguishing Agent	Carbon dioxide; foam; dry chemical
Critical Raw Materials Present in Battery in a Concentration Greater than 0.1% w/w.	Tin, Sulfuric Acid
Instructions and Safety Information	See product manual 6390-805-005
Storage and Transport Conditions	Store batteries in cool, dry, well-ventilated areas with impervious surfaces and adequate containment in the event of spills. Store under roof to protect against inclement weather. 77 Degrees F (25 degrees C)