



4X the resolution
4X the processing power

Infinite possibilities

Smart 4K features

The Stryker Connected OR Hub allows you to stream, capture, record, and print with 4K or HD.

- 4K fluorescence imaging
- 4K intelligent lighting with AutoLight
- 4K capture and record
- 4K storage
- 4K to iPad
- 4K video playback on Hub
- 4K video overlays

Secure by design: The Connected OR Hub has been updated with several features that improve network security, such as OS hardening, application whitelisting, network vulnerability scans, write protection, and data encryption.

stryker

Maximize the value of your investment

Customized financing

Through Stryker's **Flex Financial** program, we provide a range of smart alternatives designed to fit your organization's needs. Choose from numerous payment structures that can be customized to meet your budgetary needs and help build long-term financial stability.

Maintenance & support

Stryker's **ProCare** solutions are designed to decrease downtime, improve employee satisfaction, and increase revenue opportunities. We're here to answer questions, provide maintenance and support, and help your facility provide the best care possible.

Copyright © 2018 Stryker. Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker and SPY. All other trademarks are trademarks of their respective owners or holders.

Literature Number: 1000902962 Rev. A

stryker.com



Created
to connect

1688 Advanced Imaging Modalities (AIM) 4K Platform

Introducing the Connected OR Hub, the center of your smart operating room



The Connected OR Hub

Control your devices with ease

Now you have total control of the entire visualization system from the surgical field—whether that’s from the camera head or with just your voice, it’s up to you. The new intuitive interface and improved workflow matches real-world application, making control faster and easier than ever before.

Capacitive touch screen technology responds to your touch in half the time. This is the Connected OR Hub that listens and responds to your every command.

Just the way you like it

Stryker’s Surgeon Profiles allow customized settings to auto-populate before you even enter the OR, creating a consistent surgical environment.

Patient records and surgical media are available at your fingertips—you choose where and how you want to access it.

Make your best decisions

This cohesive and comprehensive collaboration between instruments and data enables a consistently outstanding surgical experience.

This means fewer variables while having access to all the tools and information required to make the best possible clinical decisions.



stryker

Our Connected OR Hub lives up to its name, connecting you to everything and everyone you need to assure a consistently outstanding **surgical experience**.



Hub to EMR

Connect to your electronic medical record with direct HL7 integration

Uni-directional and bi-directional integration



Hub to PACS

Save images to your PACS system with DICOM integration

Pre-operatively pull patient information from the modality worklist (MWL) for a more effortless experience



Hub to FTP

Wirelessly send to a secure folder on your hospital’s network



Hub to Media Asset Management

Securely edit, archive, and share your surgical images and videos



Hub to PowerShare

Access your surgical images from anywhere, at any time from the cloud



Hub to iPad

Use Stryker’s My Patient Hub iPad application to share custom post-operative reports and images with patients