



stryker

ENT

product catalog

Balloon dilation

Navigation

Bioresorbables

Cryotherapy

Nasal implants

Surgical instruments

Visualization

We are Stryker

**Together with our customers,
we are driven
to make healthcare better**

**Our mission is at the heart of what we do
and believe.**

We collaborate with our customers to develop innovative products and services that ultimately improve the lives of patients. We understand that in today's challenging healthcare environment, our customers need solutions that improve quality and efficiency in the operating room and, ultimately, the patient experience. These challenges are creating opportunities for us to make healthcare better, for our customers and the patients they serve.

Stryker is one of the world's leading medical technology companies, and, together with our customers, we are driven to make healthcare better. Stryker offers a diverse array of innovative products and services in ENT, Orthopaedics, Medical and Surgical, Neurotechnology and Spine that help improve patient outcomes. Stryker is active in over 100 countries around the world.

Follow @StrykerENT



Balloon dilation

XprESS ENT dilation system	03
Audion ET dilation system	04

Navigation

Stryker ENT navigation system and accessories	05
Stryker ENT navigation system instruments	07

Bioresorbables

NasoPore bioresorbable nasal dressing	09
HemoPore nasal dressing	10
XeroGel nasal/epistaxis pack	10
OtoPore bioresorbable ear dressing	11

Cryotherapy

ClariFix cryotherapy device	12
-----------------------------	----

Nasal implants

LATERA absorbable nasal implant system	13
--	----

Surgical instruments

ESSx microdebrider	14
MiniFESS surgical instrumentation	15
Reinforced anesthesia needle	16
Suction and irrigation tools	17


Visualization

FocESS image solutions	18
FocESS sinusscopes	19

Balloon dilation

XprESS[®]

ENT dilation system



Product	Ordering no.	Description
XprESS LoProfile system	LPLF-106	6 x 20 mm, 1 system
XprESS LoProfile system	LPLF-107	7 x 20 mm, 1 system
XprESS LoProfile system	LPLF-205	5 x 8 mm, 1 system
XprESS LoProfile system	LPLF-206	6 x 8 mm, 1 system
XprESS Ultra system	ULF-106	6 x 20 mm, 1 system
XprESS Ultra system	ULF-206	6 x 8 mm, 1 system

1 system includes: 1 balloon dilation device, 1 inflation syringe, 1 extension line, 1 bending tool, 1 PathAssist LED Light Fiber (5110-001)

HCPCS C1726 Catheter, balloon dilation, non-vascular

Balloon dilation

Audion[™]

ET dilation system




Product	Ordering no.	Description
Audion ET dilation system	AET-100	1 system

1 system includes: 1 balloon dilation device, 1 inflation syringe, 1 extension line, 1 inflation lock


HCPCS C1726 Catheter, balloon dilation, non-vascular

Stryker ENT navigation system and accessories

Scopis ENT software



Scopis ENT software with TGS (Target Guided Surgery)



Electromagnetic navigation unit

Product	Ordering no.	Description
Stryker ENT navigation system		
Electromagnetic navigation unit	8000-010-003	High-performance navigation PC with electromagnetic measuring technology
Scopis ENT software	8000-021-001	Basic software for navigation in ENT procedures
Scopis ENT software with TGS	8000-021-002	Premium target guided software with advanced planning features, virtual endoscopy and augmented reality technology
Scopis planning software	8000-021-005	Planning software for preoperative surgical planning from a laptop or workstation

Stryker ENT navigation system and accessories (cont.)




Field generator



Field generator mounting arm

Product	Ordering no.	Description
Navigation platform accessories		
Field generator	8000-010-004	Used for the electromagnetic position measurement of patient tracker and navigated instruments
Field generator mounting arm	8000-010-005	Flexible mounting arm to position the field generator
Universal mounting arm	8000-010-008	Smaller arm for flexible positioning of the field generator
Headrest - adapter A package	8000-011-005	TGS universal headrest, mounting arm, adapter A, 0.25" x 2.0"
Headrest - adapter B package	8000-012-005	TGS universal headrest, mounting arm, adapter B, 0.25" x 2.5"
Headrest - adapter C package	8000-013-005	TGS universal headrest, mounting arm, adapter C, 0.25" x 3.25"
Field generator stand	8000-014-005	Stand with mounting arm to position the field generator
Equipment cart pro	8000-030-002	Includes: four shelves and a storage drawer
Medical keyboard	8000-030-010	U.S./international keyboard, wired, number pad, USB
Mouse, wired	8000-030-020	Wired, USB, three buttons, scroll function
Mouse, wireless	8000-030-021	Wireless, USB, scroll function

Stryker ENT navigation system instruments



TGS guidewire

Pointers and suction tubes

Instrument clamp electromagnetic sphere

Instrument clamp electromagnetic universal

Product	Ordering no.	Description
Navigation instruments, single use		
TGS guidewire	8000-060-009	A navigation instrument compatible with XprESS LoProfile ENT dilation system
Navigation instruments, maximum 10 uses		
Patient tracker electromagnetic	8000-040-002	Used to locate the patient within the EM field
Precision pointer electromagnetic	8000-050-001	1.5 mm, malleable
Suction tube, Frazier electromagnetic	8000-050-005	3 mm, straight suction
Suction tube, Eicken electromagnetic	8000-050-006	3 mm, curved suction
Endoscope tracker electromagnetic	8000-060-001	Navigation suitable for 0°, 30° and 45° endoscopes with 4 mm diameter and 175 mm length
Calibration body electromagnetic	8000-060-002	Calibration of suitable conventional surgical instruments
Universal tracker electromagnetic	8000-060-006	Allows navigation of rigid conventional surgical instruments in combination with an instrument clamp
Instrument clamp electromagnetic sphere	8000-060-020	Adapter to track rigid conventional surgical instruments between 2-4.6 mm
Instrument clamp electromagnetic universal	8000-060-021	Universal adapter to track rigid conventional surgical instruments, 3-6 mm

Stryker ENT navigation system instruments (cont.)



Instrument clamp, forceps

Instrument clamp, 2-6 mm

Instrument clamp, 6-10 mm

Instrument clamp, 10-16 mm

Product	Ordering no.	Description
Navigation instruments		
Patient tracker electromagnetic	8000-040-001	Disinfectable, used to locate the patient within the EM field
Registration pointer electromagnetic	8000-050-003	Disinfectable, used for surface registration
Navigation tool extension cable	8000-050-011	Extension cable for navigation tools
Instrument clamp, forceps	8000-060-010	For attaching the universal tracker to rigid conventional surgical instruments
Instrument clamp, 2-6 mm	8000-060-011	For attaching the universal tracker to rigid conventional surgical instruments between 2-6 mm
Instrument clamp, 6-10 mm	8000-060-012	For attaching the universal tracker to rigid conventional surgical instruments between 6-10 mm
Instrument clamp, 10-16 mm	8000-060-013	For attaching the universal tracker to rigid conventional surgical instruments between 10-16 mm
Patient tracker tabs	8000-100-001	Adhesive pads for fixation of patient trackers electromagnetic, 100 per unit

NasoPore[®]

bioresorbable nasal dressings



Product	Ordering no.	Description
NasoPore dressing	5400-010-004	Standard, 4 cm, 8-pack
NasoPore dressing	5400-010-008	Standard, 8 cm, 8-pack
NasoPore dressing	5400-020-004	Firm, 4 cm, 8-pack
NasoPore dressing	5400-020-008	Firm, 8 cm, 8-pack
NasoPore FD dressing	5400-020-108	Fast degrading, 8 cm, 8-pack
NasoPore dressing	5400-030-004	Extra firm, 4 cm, 8-pack
NasoPore dressing	5400-030-008	Extra firm, 8 cm, 8-pack
NasoPore dressing	5400-212-004S	Standard, 4 cm, 2-pack
NasoPore dressing	5400-222-004S	Firm, 4 cm, 2-pack

HCPCS C1765

Cost of bleeding control/adhesion prevention/insertion of packing material, etc., is already incidental to the endoscopic sinus surgery procedures.
C1765 should be reported for XeroGel/NasoPore dressing when used with BSD or FESS procedures

HemoPore[®] and XeroGel[®]

hemostatic dressing with chitosan



Product	Ordering no.	Description
Hemopore nasal dressing		
HemoPore dressing	5400-020-208	8 cm, 8-pack
HemoPore dressing	5400-020-208S	8 cm, 2-pack
XeroGel nasal/epistaxis pack		
XeroGel pack	XG-102	4.0 L x 2.4 W x 0.3 cm thick, box of 2 individual sterile units
XeroGel pack	XG-108	4.0 L x 2.4 W x 0.3 cm thick, box of 8 individual sterile units

HCPCS C1765

Cost of bleeding control/adhesion prevention/insertion of packing material, etc., is already incidental to the endoscopic sinus surgery procedures.
C1765 should be reported for XeroGel/NasoPore dressing when used with BSD or FESS procedures

Bioresorbable dressings

OtoPore®

bioresorbable ear dressing




Product	Ordering no.	Description
OtoPore dressing	5400-010-000	Cylinder, standard, 2 x 1 cm, 8-pack
OtoPore dressing	5400-020-000	Cylinder, firm, 2 x 1 cm, 8-pack
OtoPore dressing	5400-020-100	Square, firm, 2.0 L x 4.0 W x 0.15 cm thick, 8-pack

Cryotherapy

ClariFix®

cryotherapy device



Product	Ordering no.	Description
ClariFix system	CFX-2000	1 system

1 system includes: 1 ClariFix cryotherapy device, 2 cryogen canisters and 2 sterile canister caps

HCPCS C9771 Nasal/sinus endoscopy, cryoablation nasal tissue(s) and/or nerve(s), unilateral or bilateral

LATERA[®]

absorbable nasal implant system



LATERA delivery device



LATERA implant tray and absorbable nasal implants



Product image is representative only and will include size markings on both the delivery device and implant tray.



Alar retractor

Product	Ordering no.	Description
LATERA implant system	LATSYS20	20 mm, 1 system
LATERA implant system	LATSYS24	24 mm, 1 system
Alar retractor	AR-210	Used to control and retract the ala during the implantation procedure

1 system includes: 1 implant delivery device, 2 absorbable nasal implants

HCPCS C1889 Implantable/insertable device for a device intensive procedure, not other classified

ESSx[®]

microdebrider



ESSx microdebrider

Product	Ordering no.	Description
Microdebrider individual components		
Irrigation cassettes	5290-075-000	10-pack
ESSx microdebrider	5290-601-200	Handpiece
Microdebrider blades		
Aggressive blade	5290-011-000	4.0 mm, 11 cm, 5-pack
Sharp edge blade	5290-340-000	4.0 mm, 8 cm, 5-pack
Aggressive serrated blade, 40°	5290-435-100	3.5 mm, 11 cm, 5-pack
Aggressive serrated blade, 40°	5290-440-100	4.0 mm, 11 cm, 5-pack
Jaguar blade	5290-530-000	3.5 mm, 8 cm 5-pack
Jaguar blade	5290-540-000	4.0 mm, 8 cm, 5-pack
Aggressive serrated blade	5290-628-000	2.5 mm, 11 cm, 5-pack
Aggressive serrated blade	5290-635-000	3.5 mm, 8 cm, 5-pack
Aggressive serrated blade, 60°	5290-635-100	3.5 mm, 8 cm, 5-pack
Aggressive blade	5290-638-000	3.5 mm, 8 cm, 5-pack
Aggressive serrated blade, 60°	5290-640-100	4.0 mm, 11 cm, 5-pack
Aggressive serrated blade	5290-645-000	4.0 mm, 8 cm, 5-pack
Hooded rounded bur	5290-647-000	3.0 mm, 8 cm, 5-pack
Aggressive blade	5290-648-000	4.0 mm, 8 cm, 5-pack
Angled aggressive blade	5290-740-100	4.0 mm, 11 cm, 5-pack
Convex angled aggressive blade	5290-740-200	4.0 mm, 11 cm, 5-pack
Angled aggressive blade, 12°	5290-744-100	4.0 mm, 11 cm, 5-pack
Angled aggressive blade, 60°	5290-760-100	4.0 mm, 11 cm, 5-pack
Accessories		
NSE footswitch	5400-007-000	Four feature buttons
CORE 2 console	5400-052-000	Single control center for a multitude of powered instruments
CORE footswitch	5402-007-000	Raised toe loop 270° open access


MiniFESS®

surgical instrumentation

Product	Ordering no.	Description
MiniFESS turbinate forceps	CS-101	110 mm working length, 5 x 19 mm jaws
MiniFESS Blakesley thru-cut forceps	CS-102	110 mm working length, straight, 2 x 5 mm jaws
MiniFESS Blakesley thru-cut forceps	CS-103	110 mm working length, 45°, 2 x 5 mm jaws
MiniFESS Takahashi forceps	CS-104	Straight, 2.5 x 8 mm jaws
MiniFESS sickle knife	CS-105	Curved, 10 mm
MiniFESS Blakesley forceps	CS-106	Straight, 110 mm working length, 2.5 x 7 mm jaws
MiniFESS Blakesley forceps	CS-107	45°, 110 mm working length, 2.5 x 7 mm jaws
MiniFESS Blakesley forceps	CS-108	90°, 110 mm working length, 2.5 x 7 mm jaws
MiniFESS sphenoid seeker/freer	FS-100	Sphenoid seeker/freer
MiniFESS maxillary seeker	MS-120	120° angle
MiniFESS maxillary seeker	MS-135	135° angle
MiniFESS maxillary seeker set	MS-255	Includes: both 120° angle and 135° angle
MiniFESS Light Seeker	IS-100	78°, 1.5 mm ball tip, includes standard ACMI connection plus 1 Wolf and 1 Storz connector adapter
Portable red LED light source	PRLS-100	Portable, non-sterile, reusable light source to be used with the MiniFESS Light Seeker, includes 2 batteries and 1 battery charger
MiniFESS surgical instrument set	CST-1000	Includes: 3 Blakesley forceps (straight, 45°, 90°), 2 Blakesley thru-cut forceps (straight, 45°), 1 Takahashi forceps, 1 turbinate forceps, 1 sickle knife, and 1 sterilization tray
MiniFESS seeker set	PASS-100	Includes: 2 maxillary seekers (120° and 135° angle), 1 sphenoid seeker/freer and 1 MiniFESS Light Seeker
MiniFESS suction set	STK-00	Includes: MiniFESS suctions: straight, 70° and 90°
MiniFESS suction, straight	ST-0	3 mm, straight
MiniFESS suction, 70°	ST-70	3 mm, 70°, curved
MiniFESS suction, 90°	ST-90	3 mm, 90°, curved
MiniFESS sterilization tray	TRA-101	Autoclavable; holds 8 finger-grip and 6 seeker-style instruments

RAN

(reinforced anesthesia needle)



Product	Ordering no.	Description
Reinforced anesthesia needle	RAN-027-2	27 G, 3.5", 2-pack (only sold in a bundle configuration)
Reinforced anesthesia needle	RAN-027-5	27 G, 3.5", 5-pack

Suction and irrigation tools




Product	Ordering no.	Description
Cliq suction pump	LCP-AS-1001D	Includes: 800 mL collection canister, microbial filter, intermediate tubing and suction tubing
Cyclone system	CC-100	Includes: 1 Cyclone device, 1 syringe, 1 bending tool

FocESS[®] image solutions



Product	Ordering no.	Description
FocESS HD wireless camera	WCS-100	Includes: 1 receiver, 1 battery charger and 2 batteries
Monitor, 32"	OPM-100-32	32" HD monitor


FocESS[®] sinusscopes



Product	Ordering no.	Description
FocESS sinuscope, 0°	ES-0	3 mm, 0°, autoclavable
FocESS sinuscope, 30°	ES-30	3 mm, 30°, autoclavable
FocESS sinuscope, 45°	ES-45	3 mm, 45°, autoclavable
FocESS sinuscope, 70°	ES-70	3 mm, 70°, autoclavable
FocESS sterilization tray	TRA-102	Sterilization tray for FocESS sinusscopes; autoclavable; holds 3 sinusscopes and 6 light post adapters

Unmatched support

Stryker provides the technology, training and service you need to help you improve your patient outcomes. With Stryker, you gain access to a comprehensive suite of solutions to help make your vision of patient care a reality. From technology to training, from reimbursement tools to patient education, Stryker is there to help.



Reimbursement

Today's healthcare environment is ever-evolving, and keeping up is critical for your practice. That's why we offer a variety of services and support tools to help you navigate the intricacies of reimbursement.

Contact our reimbursement hotline

Talk live with an expert to address your questions.

P 844 842 4232 | ENT-reimbursement@stryker.com



Medical Education

Whether you are in the office, the OR, or at home, our goal is to provide a well-rounded educational experience. We are committed to providing superior training on your schedule through diverse courses led by distinguished faculty in both live and virtual formats.


For upcoming events, contact our Medical Education team at ENTMedEd@stryker.com.



Flex Financial

Our payment options deliver smart alternatives designed to fit your business needs. We offer flexibility beyond a cash purchase—enabling you to secure the solutions you need for your practice, patients and staff. No other company provides strategic financial consulting combined with our breadth of products, meeting both your financial and equipment needs.

Contact your Sales Representative for more information about our flexible payment options.



ProCare Services

Safety and reliability are two of the most important reasons to use an OEM. As your trusted partner, we offer you and your new capital our expert and highly trained technicians, proprietary product knowledge, loaner equipment to maximize uptime, and the flexible and unique service solutions you need.

Contact your Stryker Sales Rep today to get a customized ProCare Services plan for all of your Stryker equipment.

ENT ordering information

To order, contact Customer Service | P 800 253 3210 | F 877 790 7772 | inst.stryker.cs@stryker.com

ENT

See Instructions for Use - visit: stryker.com/us/en/ent/ifu.html

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Audion, Building Blocks, ClariFix, Cyclone, ESSx, Flex Fincial, FocESS, HemoPore, LATERA, Light Seeker, MiniFESS, NasoPore, OtoPore, PathAssist, ProCare Services, Scopis, Stryker, TGS, XeroGel and XprESS. All other trademarks are trademarks of their respective owners or holders.

ENT-GSNPS-PRDC-655033 | 07/2023 | Copyright © 2023 Stryker

Stryker ENT
3600 Holly Lane North,
Suite 40
Plymouth, MN 55447
O 866 620 7615
F 763 463 1599

stryker.com