


# EasyStep™ Staple

## Simplifying the stepped staple experience

EasyStep is a superelastic nickel-titanium (Ni-Ti) alloy staple<sup>1</sup> with new instrumentation that is designed to simplify your stepped staple experience. The staple requires no heat activation<sup>1</sup> to provide compression between the barbed legs. The stepped staple is intended for bone fragment and osteotomy fixation of the foot.<sup>1</sup>



# EasyStep

## Stepped staple system

### Simplifying the stepped staple experience

#### Range

Offering five step heights to allow various step depths

#### Simple

Single incision procedure: prepare osteotomy and implant through the same incision

#### Flexible

Adjustable drill guide and forceps ensure the correct size implants is always used


#### Efficient

No heat activation required to provide compression between barbed legs

#### Adjustable forceps



#### Adjustable drill guide



Visit [footankle.stryker.com](http://footankle.stryker.com) or email us at [footandankle@stryker.com](mailto:footandankle@stryker.com).

#### References

1. EasyStep Op Tech ES-ST-2, 01-2016.
2. EasyStep System Components ES-SCS-2,01-2016.
3. Assembly/disassembly guide for adjustable forceps ES-ADI-1,06-2016.
4. Assembly/disassembly guide for adjustable drill guide ES-ADI-2,06-2016.

## Foot & Ankle

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: EasyStep, Stryker. All other trademarks are trademarks of their respective owners or holders.

**Content ID: ES-SS-1, 08-2016**

Copyright © 2016 Stryker



#### Manufacturer:

Stryker GmbH  
Bohnackerweg 1  
2545 Selzach  
Switzerland

[stryker.com](http://stryker.com)