


Lesser Metatarsal Head Implant

System features

- Thin, low-profile design for minimal bone resection
- Stem offset dorsally for anatomically correct alignment in medullary canal
- Rectangular shape base minimizes chance of rotation
- Vertical cut angle reduces potential for dorsal impingement

Implant features

- Single-piece, sterile-packed implant
- Highly polished cobalt chrome articular surface
- Plasma spray titanium coating on stem



Part #	Alt. part #	Description	Head height A	Head width B	Stem length C
62101001	LMH-0001	Implant size 1	10.5mm	8mm	12mm
62101002	LMH-0002	Implant size 2	13mm	10mm	12mm
62101003	LMH-0003	Implant size 3	13mm	12mm	12mm