

# PROstep<sup>®</sup> MIS Lapidus

Powerful fixation stability through a minimal incision




# PROstep<sup>®</sup> MIS Lapidus

The PROstep MIS Lapidus system is an internal fixation system intended for minimally invasive reduction of hallux valgus deformity and subsequent fusion of the first metatarsal cuneiform joint.

## MIS Joint Preparation


TMT Joint Distractor-Compressor controls access to joint for prep and powerful compression.

Arthrodesis burrs for MIS joint prep designed to preferentially remove cartilage to maintain bone length.



## Triplanar Reduction

IM angle reduction and de-rotation of the metatarsal, including plantar dorsal control. Features an ergonomic distal design.



## Fixation Stability

Three screw construct positioning designed to create a tension band effect for biomechanical stability.



## Ordering information

### PROstep MIS Lapidus Instrument (57LPKIT1)

Product no.	Description
57LPREDC	Reduction Clamp
57LPDT20	T20 Driver, Reusable
57LPDCPG	Dist-Comp Placement Guide
57LPDSCP	Distractor - Compressor
57LPDCSL	Dist - Comp Sleeve
57LPSTGL	Targeting Guide - L
57LPSTGR	Targeting Guide - R
57LPTP50	Chamfered Tissue Protector
57LPKTGR	Keystone Targeting Guide - R
57LPKTGL	Keystone Targeting Guide - L
58871010	Ratcheting Handle
57LPTRAY	Tray

### PROstep MIS Lapidus Consumables Kit (57LPKITA)

Product no.	Description
57S1MI07	PROstep Inst Pack
57KWPACK	Temp Wire Pack
58PM2SLV	Irrigation Sleeve for PM2
58TPXSLV	Irrigation Sleeve for TPX
58CC3012	Burr Cartilage Cyl 3x12
58SC3012	Burr Sculpting Cyl 3x12
58CF4008	Burr Cartilage Flame 4x8
58RSPH40	Burr Cortical Sphere 4mm
58RW3113	Burr Cortical Wedge 3.1x13
58FENAWL	Fenestration Awl
57DRDEP5	Drill and Depthing Pack 5mm
DSDS1014S	1.4mm MICA Wire
57S00030	Drill for 4mm MICA
57S02025	Driver for 4mm MICA
57S05050	5x50mm Chamfer Screw
57S05052	5x52mm Chamfer Screw
57S05054	5x54mm Chamfer Screw
57S05056	5x56mm Chamfer Screw
57S05058	5x58mm Chamfer Screw

### PROstep MIS Lapidus Consumables Kit (57LPKITA) continued

57S05060	5x60mm Chamfer Screw
57S05065	5x65mm Chamfer Screw
57S05070	5x70mm Chamfer Screw
57S05075	5x75mm Chamfer Screw
57S05080	5x80mm Chamfer Screw
57S34020	4x20mm MICA Screw
57S34042	4x42mm MICA Screw
57S34022	4x22mm MICA Screw
57S34044	4x44mm MICA Screw
57S34024	4x24mm MICA Screw
57S34046	4x46mm MICA Screw
57S34026	4x26mm MICA Screw
57S34048	4x48mm MICA Screw
57S34028	4x28mm MICA Screw
57S34050	4x50mm MICA Screw
57S34030	4x30mm MICA Screw
57S34052	4x52mm MICA Screw
57S34032	4x32mm MICA Screw
57S34054	4x54mm MICA Screw
57S34034	4x34mm MICA Screw
57S34056	4x56mm MICA Screw
57S34036	4x36mm MICA Screw
57S34058	4x58mm MICA Screw
57S34038	4x38mm MICA Screw
57S34060	4x60mm MICA Screw
57S34040	4x40mm MICA Screw
57S34030	4x30mm MICA Screw
57S34032	4x32mm MICA Screw
57S34034	4x34mm MICA Screw
57S34036	4x36mm MICA Screw
57S34038	4x38mm MICA Screw
57S34040	4x40mm MICA Screw
57S34042	4x42mm MICA Screw
57S34044	4x44mm MICA Screw
57S34046	4x46mm MICA Screw
57S34048	4x48mm MICA Screw
57S34050	4x50mm MICA Screw

## Foot & Ankle

Visit us at [footankle.stryker.com](https://footankle.stryker.com)

The information contained in this document is intended for healthcare professionals only. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: PROstep, Stryker. All other trademarks are trademarks of their respective owners or holders.

AP-016877, 08-2023

Copyright © 2023 Stryker



**Manufacturer:**  
Wright Medical Technology  
1023 Cherry Road  
Memphis, TN 38117  
800 238 7117  
901 867 9971

[www.wright.com](https://www.wright.com)