

# VariAx<sup>®</sup> 2

## Forefoot fracture plates



**Plating** that fits.

# VariAx 2

## Forefoot fracture plates

### Fits for patients

- Designed for an enhanced fit
- Low profile design/1.5mm thickness

### Fit for procedures


- Variable angle locking (15°)
- Type II anodized titanium alloy
- Versatility of plate choices

### Fit for you

- Simple instrumentation
- Flexible tray configuration




Slim/T8 plates



2.4mm


2.7mm



14 – 49mm

2 hole – 7 hole

Thickness: 1.5mm



30 – 64mm

2 hole – 7 hole

Thickness: 1.5mm



Broad/T10 plates



2.7mm


3.5mm



19 – 60mm

2 hole – 7 hole

Thickness: 1.5mm



36 – 77mm

2 hole – 7 hole

Thickness: 1.5mm

Visit [footankle.stryker.com](http://footankle.stryker.com) or email us at [footandankle@stryker.com](mailto:footandankle@stryker.com).

## Foot & Ankle

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, VariAx. All other trademarks are trademarks of their respective owners or holders.

This document is only approved for use in the USA.

Content ID: VAX-SS-10, 06-2017  
Copyright © 2017 Stryker

**Manufacturer:**  
Stryker GmbH  
Bohnackerweg 1  
2545 Selzach,  
Switzerland

[stryker.com](http://stryker.com)