

VariAx[®] 2

MIS Calcaneus

Continuing to advance our plating portfolio



VariAx[®] 2 MIS Calcaneus

Plates

VariAx 2 MIS Calcaneal Plates

MIS Calcaneal Plate,
U5, 5 holes



MIS Calcaneal Plate,
U6, 6 holes



MIS Calcaneal Plate,
S7, 7 holes



MIS Calcaneal Plate,
M8, 8 holes



MIS Calcaneal Plate,
L9, 9 holes



Screws

The plates within the VariAx 2 MIS Calcaneus are compatible with VariAx 2 screws.

Locking Non-locking Locking Non-locking



**Ø3.5mm screws
T10 interface**



**Ø2.7mm screws
T10 interface**



Key features

- 5 plating options
- Non side specific – interchangeable between the left and right calcanei
- VariAx 2 T10 instrumentation for Ø2.7mm/3.5mm screws
- Elevator, soft tissue instrument doubles as an x-ray template
- Targeting guide for MIS approach



Foot & Ankle

Visit us at footankle.stryker.com

The information contained in this document is intended for healthcare professionals only. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: VariAx, Stryker. All other trademarks are trademarks of their respective owners or holders. This document is only approved for use in the USA.

VAX-SS-32, 09-2023
Copyright © 2023 Stryker

Manufacturer:
Stryker GmbH
Bohnackerweg 1
2545 Selzach,
Switzerland
stryker.com