

**stryker**

Signature Portfolio  
Elite Attachments

**Confidence**  
extended



Longer, concave design for improved line of sight and fine-edge knurling for better grip and more options for positioning your hand

Snap-Connect for easy and consistent assembly/disassembly


Refined attachment design for improved balance and maneuverability

Easily identifiable marks help surgeon and staff adjust exposure to desired level

Patented multi-notch telescoping cutting tools enable convenient customization of exposure setting for specific user preferences and applications

- Adjustable up to five exposure levels (up to 10mm)
- Compatible with any Elite attachment
- 5:1 inventory reduction
- Streamlines set up and packaging waste
- Reduces opening multiple cutting accessories for different attachments

Carbide round burs have patented multi-flute design



Elite Attachments

# Efficient. Extensive. Essential.

When you're in the O.R., having access to just the right tools at your fingertips is essential. Our Elite attachments put versatility and functionality in your hands with a natural balance and grip that enhances your feel and maneuverability. Paired with our extensive selection of industry-unique telescoping cutting accessories, the Elite line gives you the ability to expand your reach by adjusting up to five exposure levels without the need to change your attachment or bur. Streamlining both ease of use and your inventory items of up to 5:1, Elite delivers efficiency all around.

## New! Irrigation Sleeves

- Deliver precise irrigation by simply sliding over your 7cm or 12-20cm Elite attachment for a complete solution
- Low profile with ergonomic design facilitates site access and visibility
- Sterile, disposable and compatible with CORE Console disposable irrigation cassettes



5407-120-XXX

## New! Elite LX Cutting Accessories

Our recent addition of 39 new Elite Cutting Accessories expands our collection of burs to more than 300 styles and options, allowing you to customize with confidence. Among them are nine carbide accessories — including the carbide matchhead — that reduce the chatter level over steel by 30% for added precision and control.<sup>1</sup>

Represented below are the new style types.

3 precision acorn;  
dimensions 5.0 – 7.5mm



3 router;  
dimensions 1.5 – 2.3mm



2 tool steel match head;  
dimensions 1.7 – 2.2mm



2 diamond match head;  
dimensions 1.7 – 2.2mm



12 diamond;  
dimensions 0.70 – 6.0mm



10 round;  
dimensions 0.70 – 7.5mm



1 carbide match head;  
dimension 3.0mm




8 carbide round;  
dimensions 2.0 – 6.0mm





## I.D. Touch Software: Define your style

When used with either the  $\pi$ drive or  $\pi$ drive+ motors, your cutting accessories can be optimized even further using I.D. Touch Software on the CORE 2 Console instrument driver. This industry exclusive software enables you to adjust torque from 0-100% to customize how your electric drill responds when pressure is applied during use — from low to high sensitivity — redefining performance to your liking. Your NSE sales representative can guide you through the options and save the performance and feel that's right for you in user preferences.





## Signature Portfolio Elite Attachments and Cutting Accessories

To order our Elite attachments and cutting accessories – or to learn more about our Signature Portfolio – please visit [neurosurgical.stryker.com](http://neurosurgical.stryker.com) or call 800 253 3210.

Ordering information		
Part number	Product description	Color band
5407-120-450	7cm Elite Straight Attachment	Black
5407-120-452	7cm Elite Angled Attachment	Black
5407-120-470	12cm Elite Straight Attachment	Black
5407-120-472	12cm Elite Angled Attachment	Black
5407-120-480	14cm Elite Straight Attachment	Black
5407-120-482	14cm Elite Angled Attachment	Black
5407-120-484	14cm Elite Hooded Attachment	Black
5407-120-485	17cm Elite Straight Attachment	Black
5407-120-486	17cm Elite Angled Attachment	Black
5407-120-490	20cm Elite Straight Attachment	Black
5407-120-492	20cm Elite Angled Attachment	Black
5407-120-050	9cm Heavy Duty Straight Attachment	Blue
5407-120-052	9cm Heavy Duty Angled Attachment	Blue
5407-120-070	14cm Heavy Duty Straight Attachment	Green
5407-120-072	14cm Heavy Duty Angled Attachment	Green

Signature Portfolio Elite Irrigation		
Part number	Product description	
5407-010-450	Irrigation sleeves for use with 7cm (black) Elite attachments 5407-120-450      5407-120-452	
5407-010-470	Irrigation sleeves for use with 12-20cm (black) Elite attachments 5407-120-470      5407-120-472      5407-120-480 5407-120-482      5407-120-484      5407-120-485 5407-120-486      5407-120-490	
5407-010-270	7cm Elite Irrigation clip	
5407-010-272	12-20cm Elite Irrigation clip	
5407-010-200	9cm HD Irrigation clip	
5100-010-200	Universal irrigation	

1. Stryker internal data on file, Kalamazoo, Michigan, November 2016.

## Neurosurgical

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: CORE, I.D. Touch,  $\pi$ drive and Stryker. All other trademarks are trademarks of their respective owners or holders.