

stryker

The CORE FS

In step
with **you**



Thinking one step ahead

During surgery, your equipment needs to perform on demand. That's precisely why we engineered the CORE footswitch with a design that's ergonomic and easily accessible from any angle – even while standing or sitting. With the power and precision you can count on, it's always ready for your next step.



Raised toe loop:

- Enables easy positioning and movement by user
- Toe loop in stow (lowered) position allows convenient cord wrapping for storage

Non-slip rubber feet
designed to minimize sliding

Compatible with
CORE 2 Console software

Treadle design
has ideal height and spring force

Coloured light on each side:

- Shows which drill is live
- Maps to CORE 2 Console footswitch ports (blue=top port, amber=bottom port)

No raised side walls/edges
for 270° open access

Pedal raises to allow for thorough cleaning and disinfecting

Learn more

For more information, contact your NSE sales representative

Product	Reference number
CORE Footswitch	5402-007-000
CORE 2 Console (power cord included)	5400-052-000

NSE

This document is intended solely for the use of healthcare professionals.

SA healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: CORE and Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



D0000007337EN

SDL 10/2019

2019-20310

Copyright © 2019 Stryker

stryker.com



Stryker Neuro Spine ENT
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
Tel.: (+1) 269 323 7700
Fax: (+1) 800 999 3811