

**stryker**

Protection  
and power

**perfectly  
balanced**



**Signature Portfolio  
Footed Attachments**

Footed Attachments

# Putting the most advanced footed attachments in your hands

When you told us that balancing power while protecting delicate tissue is critical in neurosurgery, we listened. It's the reason we carefully engineered our Signature Portfolio Footed Attachments for the ideal combination of precision and power every time. Using your input, we refined and optimized our components to give you enhanced maneuverability and control even in the most demanding of procedures. To help improve your overall experience, we implemented our Snap-Connect color coding system to align attachments with our lineup of compatible accessories making them immediately identifiable and easy to use.



## Customize with confidence

We understand that needs change from patient to patient. That's why we offer a spectrum of Signature Portfolio routers and cutting accessories in various head types, diameters (1.1-3.0mm) and lengths (12-25mm) to create the exact functionality you want.

Match Head



Diamond Match Head



Straight Router



Tapered Router



Spiral Router



Cutting accessories lock directly into motor, allowing them to remain in place when switching attachments

New enhanced ceramic ball bearings extend attachment lifespan



Attachments and cutting accessories can be easily paired:

- During set up via color-coded attachments and cutting accessory packaging
- On the surgical field via matching part numbers

Tactile toe alignment indicator signals which way foot is pointing

Fine-edge knurling enhances grip and comfort

Alignment dots guide connection of attachments to motor


Snap-Connect for easy and consistent assembly/disassembly

Exclusive πdrive motor is highly tunable, allowing you to dial in on how you want the drill to feel



## I.D. Touch Software: Define your style

When used with either the  $\pi$ drive or  $\pi$ drive+ motors, your cutting accessories can be optimized even further using I.D. Touch Software on the CORE 2 Console instrument driver. This industry exclusive software enables you to adjust torque from 0-100% to customize how your electric drill responds when pressure is applied during use — from low to high sensitivity — redefining performance to your liking. Your Neurosurgical sales representative can guide you through the options and save the performance and feel that's right for you in user preferences.



## Signature Portfolio Footed Attachments and Cutting Accessories

To order our Footed attachments and cutting accessories – or to learn more about our Signature Portfolio – please visit [neurosurgical.stryker.com](http://neurosurgical.stryker.com) or call **800 253 3210**.

### Attachments

Part number	Product description	Color band
5407-FA1-000	12mm Footed Attachment	Green
5407-FA2-000	16mm Footed Attachment	Blue
5407-FA3-000	25mm Footed Attachment	White
5407-FA1-00R	12mm Rotating Footed Attachment	Green
5407-FA2-00R	16mm Rotating Footed Attachment	Blue
5407-FA3-00R	25mm Rotating Footed Attachment	White
5407-120-450	7cm Straight Attachment	Black
5407-NFA-8CM	8cm Non Footed Attachment	Red
5407-NFA-9CM	9cm Non Footed Attachment	White

### Cutting Accessories

5407-FA1-015	1.5mm Tapered Router	Green and red
5407-FA1-015SP	1.5mm Spiral Router	Green and red
5407-FA2-023	2.3mm Tapered Router	Blue and red
5407-FA2-023SP	2.3mm Spiral Router	Blue and red
5407-FA3-030	3.0mm Tapered Router	White
5407-FA3-030SP	3.0mm Spiral Router	White
5407-107-017	1.7mm Match Head	Red
5407-107-017D	1.7mm Diamond Match Head	Red
5407-107-022	2.2mm Match Head	Red
5407-107-022D	2.2mm Diamond Match Head	Red
5407-070-011	1.1mm Straight Router	Red
5820-070-011	1.1mm Straight Router	Black
5820-070-015	1.5mm Straight Router	Black
5820-070-020	2.0mm Straight Router	Black
5820-070-025	2.5mm Straight Router	Black
5820-070-030	3.0mm Straight Router	Black
5820-071-015	1.5mm Tapered Router	Black
5820-071-020	2.0mm Tapered Router	Black
5820-071-025	2.5mm Tapered Router	Black
5820-071-030	3.0mm Tapered Router	Black
5820-071-023	2.3mm Tapered Router	Black
5820-071-015SP	1.5mm Spiral Router	Black
5820-071-023SP	2.3mm Spiral Router	Black

NFA = Non footed attachments

## Neurosurgical

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: CORE, I.D. Touch,  $\pi$ drive and Stryker. All other trademarks are trademarks of their respective owners or holders.