

stryker



OtoPore[®]

Bioresorbable Ear Dressing

For outer and middle ear surgery

Made of bioresorbable synthetic polymer, OtoPore ear dressings provide post-op structural support in various shapes and sizes to meet the needs of different anatomical placements and user preferences. Its biocompatible formula self degrades into natural components and clears through the Eustachian tube or mucosal flow.

Predictable performance



OtoPore bioresorbable synthetic polyurethane foam is available in different densities, shapes and sizes. It aids the post-op structural support.

- Provides haemostasis through compression
- Minimises adhesions and edema
- Absorbs blood/fluid
- Separates and supports vital tissue structures
- Clears naturally
- Easy to cut and manipulate during placement
- No mixing, moistening or other prep needed

| | Ear Cylinder Standard | Ear Cylinder Firm | Ear Square Firm |
|---|------------------------------------|------------------------------------|-------------------------------------|
| Order number (QTY/8 per box) | 5400-010-000ITL (2 x 1cm) | 5400-020-000ITL (2 x 1cm) | 5400-020-100ITL (2cm x 4cm x 1.5mm) |
| Typical placement (more than one unit may be applied if needed) | Outer ear | Outer ear | Middle ear |
| Fragmentation time* | 7-15 days | 7-15 days | 14-28 days |
| Compression | 36-48 hours | 36-48 hours | 36-48 hours |
| Hydrophilic | Yes | Yes | Yes |
| Material | 100% synthetic proprietary polymer | 100% synthetic proprietary polymer | 100% synthetic proprietary polymer |

*May vary due to amount used, degree of saturation and site.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: OtoPore, Stryker. All other trademarks are trademarks of their respective owners or holders.

Literature Number: 9100-003-282 Rev. A_SSP
Copyright © 2017 Stryker

Distributed by:
Stryker Australia
8 Herbert Street,
St Leonards NSW 2065
t: +61 2 9467 1000
f: +61 2 9467 1010

Stryker New Zealand Limited
PO Box 17-136
Greenlane, Auckland 1546
T: 64 9 573 1894
F: 64 9 573 1895

www.StrykerMedEd.com