

Sonopet iQ

Ultrasonic Aspirator System

Bone cutting tips

12cm iQ Micro Claw

5500-25B-315
4/box

L 4.8in (121mm) **CL** 1.6mm **CW** 2mm



12cm iQ Claw

5500-25B-318
4/box

L 4.8in (121mm) **CL** 1.7mm **CW** 2.8mm



12cm iQ 120 Claw

5500-25B-317
4/box

L 4.8in (121mm) **CL** 2.58mm **CW** 3.15mm



12cm iQ Apex 360

5500-25B-312
4/box

L 4.8in (121mm) **CL** 2.58mm **CW** 3.15mm



20cm iQ Apex 360

5500-25B-302
4/box

L 7.9in (200mm) **CL** 2.5mm **CW** 3.15mm



18cm iQ Micro Claw

5500-25B-303
4/box

L 7.2in (182.8mm) **CL** 1.6mm **CW** 2mm



18cm iQ Apex Micro Claw

5500-25B-319
4/box

L 7.2in (182.8mm) **CL** 1.6mm **CW** 2mm



11cm iQ Knife

5500-25B-107
4/box

L 4.2in (106mm) **CL** 12.5mm **CW** 0.83mm



11cm iQ Apex Knife

5500-25B-114
4/box

L 4.2in (106mm) **CL** 12.5mm **CW** 0.83mm



L = Length **CL** = Cutting length **CW** = Cutting width

12cm iQ Micro

5500-25S-309

4/box

L 4.7in (120mm) **OD** 1.6mm **ID** 1.20mm**12cm iQ Standard**

5500-25S-307

4/box

L 4.5in (115mm) **OD** 1.92mm **ID** 1.5mm**12cm iQ Large**

5500-25S-308

4/box

L 4.7in (112.1mm) **OD** 2.49mm **ID** 2.0mm**20cm iQ Standard**

5500-25S-301

4/box

L 8.1in (205mm) **OD** 1.92mm **ID** 1.5mm**12cm iQ Barracuda**

5500-25S-310

4/box

L 4.5in (115mm) **OD** 1.92mm **ID** 1.5mm**12cm iQ Serrated Large**

5500-25S-316

4/box

L 4.6in (115.7mm) **OD** 2.51mm **ID** 2.0mm**L** = Length **OD** = Outer diameter **ID** = Inside diameter**Neurosurgical**

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Sonopet iQ and Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

D0000011657 Rev. AA
G55/PS
Copyright © 2019 Stryker
Printed in USA

Stryker
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 269 323 2742
toll free: 800 253 3210

neurosurgical.stryker.com