

Signature Portfolio Footed Attachments and Cutting Accessories

Attachments		
Part Number	Product Description	Color Band
5407-FA1-000	Footed Attachment 12mm	Green
5407-FA2-000	Footed Attachment 16mm	Blue
5407-FA3-000	Footed Attachment 25mm	White
5407-FA1-00R	Rotating Footed Attachment 12mm	Green
5407-FA2-00R	Rotating Footed Attachment 16mm	Blue
5407-FA3-00R	Rotating Footed Attachment 25mm	White
5407-NFA-8CM	Non Footed Attachment 8cm	Red
5407-NFA-9CM	Non Footed Attachment 9cm	White
Cutting Accessories		
Part Number	Product Description	Color Band
5407-FA1-015	1.5mm Tapered Router	Green & Red
5407-FA1-015SP	1.5mm Spiral Router	Green & Red
5407-FA2-023	2.3mm Tapered Router	Blue & Red
5407-FA2-023SP	2.3mm Spiral Router	Blue & Red
5407-FA3-030	3.0mm Tapered Router	White
5407-FA3-030SP	3.0mm Spiral Router	White
5407-107-017	1.7mm Match Head	Red
5407-107-017D	1.7mm Diamond Match Head	Red
5407-107-022	2.2mm Match Head	Red
5407-107-022D	2.2mm Diamond Match Head	Red
5407-070-011	1.1mm Straight Router	Red
NFA = Non footed attachments		

To order footed attachments and cutting accessories – or to learn more about our Signature Portfolio – please visit www.stryker.com/NSE or call **800-253-3210**.

NSE
Neuro Spine ENT

We're laser focused on Neuro, Spine and ENT. Our instruments, technology, education and support are all tailored to your specialties to help fulfill your unique needs and potentially make a positive difference in your day and outcomes.

- High Speed Drills (electric and pneumatic)
- Attachments and cutting accessories
- Sonopet® Ultrasonic Aspirator
- Bone Mill
- SilverGlide® Bipolar Forceps
- Microdebrider
- NasoPore® Bioresorbable Nasal Dressings
- HemoPore® Hemostatic Bioresorbable Nasal Dressing
- OtoPore® Bioresorbable Ear Dressings
- Surgeon headlights and other accessories in support of NSE specialties



Reconstructive

- Hips
- Knees
- Trauma & Extremities
- Joint Preservation
- Orthobiologics

Medical & Surgical

- Power Tools & Surgical Accessories
- Image Guided Navigation
- Endoscopy & Arthroscopy
- Integrated Communications
- Beds, Stretchers & EMS
- Sustainability Solutions

Neurotechnology & Spine

- Craniomaxillofacial
- Interventional Spine
- Neurosurgical, Spine & ENT
- Neurovascular
- Spinal Implants

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with TM designation are trademarks of Stryker. Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 9100-003-172 Rev None G55/PS

Copyright © 2015 Stryker
Printed in USA

Stryker Neuro Spine ENT
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 800 999 3811
toll free: 800 253 3210

www.stryker.com/NSE



Neuro Spine ENT



Footed Attachments




Footed



A hands-down advancement in **footed instruments**

The Signature Portfolio marks Stryker's most comprehensive high speed drill launch, and includes a new generation of footed attachments and cutting accessories. Each component has been refined and optimized based on user input to enhance maneuverability, ease-of-use and your overall experience. Here's a look at our latest footed lineup.



Portfolio **options**

To achieve the precise functionality you want, we offer a spectrum of Signature Portfolio routers and cutting accessories in various head types, diameters (1.1-3.0mm) and lengths (12-25mm). Please see the back cover for more details.

Match Head




Diamond Match Head




Straight Router



Tapered Router



Spiral Router



Cutting accessories lock directly into motor, allowing them to remain in place when switching attachments

Snap-Connect for easy and consistent assembly/disassembly

Alignment dots guide connection of attachments to motor

Attachments and cutting accessories can be easily paired:

- During set up via color-coded attachments and cutting accessory packaging
- On the surgical field via matching part numbers

Tactile toe alignment indicator signals which way foot is pointing

Fine-edge knurling enhances grip and comfort

Exclusive π drive™ motor is highly tunable, allowing you to fine tune how you want the drill to feel